

```
Attribute VB_Name = "General"
Option Explicit

'Program and DB Rev
Public Const ProgRev$ = "pp Rev 1.00"
Public Const DBRev$ = "pp MDB Template 1.04"
```

```
'Columns collection
Public PSColumns As New Collection

'Statistics XArray
Public StatArray As New XArray
```

```
'Registry keys
Public Const SKEY_PPPPS$ = "PamelaProbeDesigner"
Public Const VALUE_DB_TEMPLATE$ = "DBTemplate"
Public Const VALUE_DB_DIR$ = "DBDir"
Public Const VALUE_DB_FILE_EXT$ = "DBFileExt"
Public Const VALUE_GB_DIR$ = "GBDir"
Public Const VALUE_GB_FILE_EXT$ = "GBFileExt"
Public Const VALUE_BLAST_DIR$ = "BLASTDir"
directory
```

```
'Error codes
Enum ppErrors
    ppErrGBFileLength = 3000
    ppErrGBFileFormat = 3001
    ppErrBadRecordset = 3002
    ppErrNoProbesCreated = 3003
    ppErrCancelled = 3004
    ppErrDBFormat = 3005
    ppErrBlastVersion = 3006
    ppErrFASTAFormat = 3007
    ppErrSeqLoadDB = 3008
    ppErrDLL = 3009
End Enum
```

```
'Flags to prevent sequence and probeset selectors firing.
'Currently broken -- RowColChange events are ignored.
Public AllowDBGClick As Boolean
```

```
'Algorithm parameter sets
Public CreatePSParms As New cCreatePSParms
Public PosFilterParms As New cPosFilterParms
Public LengthFilterParms As New cLengthFilterParms
Public GCFilterParms As New cGCFilterParms
Public dGDPars As New cdGDPars
Public dGDFilterParms As New cdGDFilterParms
Public TMPars As New cTMPars
Public TMFilterParms As New cTMFilterParms
Public dGHPars As New cdGHPars
Public dGFilterParms As New cdGFilterParms
Public dGMPars As New cdGMPars
Public dGMFilterParms As New cdGMFilterParms
Public RunParms As New cRunParms
Public RunFilterParms As New cRunFilterParms
Public Progress As New cProgress
Public EntrezParms As New cEntrezParms
```

```
Public TrimParms As New cTrimParms
Public ClampFilterParms As New cClampFilterParms
Public BlastParms As New cBLASTParms
Public HomologyParms As New cHomologyParms
Public HomoFilterParms As New cHomoFilterParms
Public ArrayParms As New cArrayParms
```

```
' seqRecord Structure
' This structure holds read from one sequence file/DB/website.
' It is used when reading Genbank files (in which case all fields
' are filled in), or FASTA files (in which case only the Header and
' sequence portions are filled in).
```

```
Public Type SeqRecord
    Header As String
    Locus As String
    Accession As String
    Length As Long
    Sequence As String
End Type
```

```
Sub Main()
'Initialize thermodynamic parameters.
InitThermoParms
```

```
'Start it up!
frmMain.Show

End Sub
```

```
Attribute VB_Name = "Blast"
Option Explicit
```

```
Private Declare Function InitPH Lib "homologydll.dll"
    (ByVal DB$, ByVal NumProbes$, ByVal W$, ByVal MaxProbeLength$) As Long
```

```
Private Declare Function AddProbe Lib "homologydll.dll"
    (ByVal Probe$, ByVal ProbeLength$) As Long
```

```
Private Declare Function CalcOneHomology Lib "homologydll.dll" (ByVal
    SeqLength$) As Long
```

```
Private Declare Function SkipSequence Lib "homologydll.dll" (ByVal SeqLength$)
    As Long
```

```
Private Declare Function GetHomology Lib "homologydll.dll" (ByRef Homologys$) As
    Long
```

```
Private Declare Function TerminatePH Lib "homologydll.dll" () As Long
```

```
Public Const MaxBLASTHeaders$ = 100
```

```
Public Sub CalcHomology(Seqs(), HP As CHomologyPairs, Homology$( ), RSblast As
Recordset, PNames$)
'-----
'Function
'Calculate the homology of an array of probes to a database of sequences.
'Arguments
'Seq: The sequences.
'HP: The homology parameters.
'Homology: The homologies.
'RSblast: The table of BLAST-located homologies.
'PName: The name of the probset we are processing.
'-----
On Error Goto E
```

```
Dim NumSeqs
Dim MaxProbelen
Dim P$, S$
Dim DBFile$
Dim SR As SeqRecord
Dim LineText$
Dim NumHeaders$
Dim DoHomology As Boolean
Dim foo$
Dim N$, NumNoCheck$
Dim NoCheck$()
Dim HomologyBaseName$

'Determine the number of probes we are calculating homology for.
NumSeqs = UBound(Seq) + 1

'Find the longest probe.
MaxProbelen = 0
For P = 0 To NumSeqs - 1
    If Len(Seq(P)) > MaxProbelen Then MaxProbelen = Len(Seq(P))
Next P

'Check for the existence of the homologizer files.
HomologyBaseName = Left(HP.DB, Len(HP.DB) - 4)
foo = FileLen(HP.Path & "\" & HomologyBaseName & ".headers")
foo = FileLen(HP.Path & "\" & HomologyBaseName & ".seqs")

'Initialize the homologizer.
If (InitPH(HP.Path & "\" & HomologyBaseName & ".seqs", NumSeqs, HP.Speed,
MaxProbelen) <> 0) Then
    Err.Raise pperrdll, "Error in InitPH DLL Call."
End If

'Load the probes into the homologizer index.
For P = 0 To NumSeqs - 1
    If (AddProbe(Seq(P), Len(Seq(P))) <> 0) Then
        Err.Raise pperrdll, "Error in AddProbe DLL Call."
    End If
Next P

'Open the specified DB.
DBFile = FreeFile
Open HP.Path & "\" & HomologyBaseName & ".headers" For Input As #DBFile

'Read the number of entries expected. The progressbar for this application
'is posted with a maximum of 100, since we don't know the number of sequences
```

```
'in advance of opening the DB.
Line Input #DBFile, LineText
NumHeaders = CLng(LineText)

'Create an array of all accessions that should not be searched against.
NumNoCheck = 0
If IsGoodRS(RSblast) Then
    RSblast.MoveLast
    NumNoCheck = RSblast.RecordCount
End If
ReDim NoCheck(0 To NumNoCheck)
N = 0
Do While Not RSblast.EOF
    If Not RSblast("Use") Then
        NoCheck(N) = RSblast("Accession")
        N = N + 1
    End If
    RSblast.MoveNext
Loop
NumNoCheck = N

'Homologize all the sequences in the DB.
S = 0
Do While Not EOF(DBFile)
    S = S + 1
    If (Fix(100 * S / NumHeaders) - Progress.StopAt > 0) Then
        Progress.CheckProgress Fix(100 * S / NumHeaders)
        'Read the next DB sequence.
        SR = ReadFASTAHeader(DBFile)
        'Check the header against the bad list.
        DoHomology = True
        For N = 0 To NumNoCheck - 1
            If SR.Accession = NoCheck(N) Then DoHomology = False
        Next N
        'Do it.
        If DoHomology Then
            If (CalcOneHomology(SR.Length) <> 0) Then
                Err.Raise pperrdll, "Error in CalcOneHomology DLL Call."
            Else
                If (SkipSequence(SR.Length) <> 0) Then
                    Err.Raise pperrdll, "Error in SkipSequence DLL Call."
                End If
            End If
        End If
    End If
    'Retrieve the answers, and clean up.
    If (GetHomology(Homology(0)) <> 0) Then
        Err.Raise pperrdll, "Error in GetHomology DLL Call."
    If (TerminatePH() <> 0) Then
        Err.Raise pperrdll, "Error in TerminatePH DLL Call."
    End Sub
Exit Sub
E: Debug.Print "Error in CalcHomology"
Err.Raise Err.Number, Err.Description
End Sub
```

```

Public Function BLASTPS$(PNames$, ByVal Sequence$, BP As cBLASTParms,
MatchAccession$(), MatchScores$(), MatchExpects$())
-----
'Function
' Run BLAST on the sequence of a probeset, saving headers and scores of
' homologues.
'Arguments
' PNames: The probeset name, used to name temporary files.
' Sequence: The base sequence of the probeset.
' BP: The BLAST search parameters.
' MatchAccession: Matching accession numbers.
' MatchHeader: The full matching header.
' MatchScores: The match scores.
' MatchExpects: The expected number of matches with such scores.
'Notes
' This is a very fragile function, since it depends on the format of the
' BLAST output file (checked only for version number), and relies on every
' header containing an accession number (though one will be manufactured
' if none is found).
-----
On Error Goto E

Dim BIFilename$, BOFilename$
Dim BIFile$, BOFile$
Dim BLASTBaseName$
Dim cmd$
Dim ProcID
Dim foo
Dim LineText$
Dim H$
Dim CBWhere$
Match Header.

'Check for the existence of all the BLAST files. If files are absent, an error
will occur.
BLASTBaseName = Left(BP.DB, Len(BP.DB) - 4)
foo = FileLen(BP.Path & "\blastall.exe")
foo = FileLen(BP.Path & "\ " & BLASTBaseName & ".nin")
foo = FileLen(BP.Path & "\ " & BLASTBaseName & ".nhr")
foo = FileLen(BP.Path & "\ " & BLASTBaseName & ".nsq")

'BLAST file names.
BIFilename = PNames & ".fsa"
BOFilename = PNames & ".out"

'Open the BLAST input file.
BIFile = FreeFile
Open BP.Path & "\ " & BIFilename For Output As BIFile

'Write out the header.
Print #BIFile, "> PNames = " & PNames

'Write out the sequence.
Do While Len(sequence) > 60
Print #BIFile, Mid$(sequence, 1, 60)
Sequence = Right(sequence, Len(sequence) - 60)
Loop

```

```

Print #BIFile, Sequence
'Close input file.
Close BIFile

'Put together the command string to execute.
BOFilename = PNames & ".out"
cmd = "cmd /c cd " & BP.Path
cmd = cmd & " & blastall -p blastn -i " & BIFilename & " -d " & BLASTBaseName &
"-o " & BOFilename

'Run it and wait.
ProcID = Shell(cmd, vbHide)
foo = WaitOnProgram(ProcID)

'Open the output file.
BOFile = FreeFile
Open BP.Path & "\ " & BOFilename For Input As #BOFile

'Check for correct version
Line Input #BOFile, LineText
If Left(LineText, 12) <> "BLASTN 2.0.2" Then Err.Raise ppErrBlastVersion

'Move down to MatchExpects.
On Err Goto ErrNoSignificantMatch
Do
Line Input #BOFile, LineText
Loop Until Instr(LineText, "sequences producing significant alignments:") <> 0
On Error Goto E

'Strip one more line
Line Input #BOFile, LineText

'Read match information, until a blank line.
Line Input #BOFile, LineText
H = 0
Do
LineText = Right(LineText, Len(LineText) - 68)
MatchScores(H) = Cdbl(Left(LineText, 4))
LineText = Right(LineText, Len(LineText) - 4)
If IsNumeric(LineText) Then
MatchExpects(H) = Cdbl(LineText)
Else
MatchExpects(H) = 0#
End If
H = H + 1
Line Input #BOFile, LineText
Loop Until Len(LineText) = 0

'We will return the number of matches.
BLASTPS = H - 1

'Read as many headers as we found matches.
For H = 0 To BLASTPS - 1
Do
'Strip lines until we see ">".
Line Input #BOFile, LineText

```

```
Loop Until Left(LineText, 1) = ">"
```

```
'Push this line onto the header.  
MatchHeader(H) = LineText
```

```
'keep pushing on left-justified lines until we see "Length".
```

```
Do  
    Line Input #BOFile, LineText  
    LineText = LTrim(LineText)  
    If (Left(LineText, 6) = "Length") Then Exit Do  
    MatchHeader(H) = MatchHeader(H) & " " & LineText  
Loop Until False
```

```
'Find the accession number.  
GBWhere = Instr(MatchHeader(H), "/"gb=")
```

```
If (GBWhere = 0) Then  
    MatchAccession(H) = "200000"  
Else  
    GBWhere = GBWhere + 4  
Do
```

```
1)    MatchAccession(H) = MatchAccession(H) & Mid(MatchHeader(H), GBWhere,  
        GBWhere = GBWhere + 1  
    Loop Until Mid(MatchHeader(H), GBWhere, 1) = " "  
    End If  
Next H
```

```
'Remove the BLAST files.
```

```
Close #BOFile  
Kill BP.Path & "\" & BFileName  
Kill BP.Path & "\" & BOFileName  
Exit Function
```

```
ErrNoSignificantMatch:
```

```
'No match was found, so we've skipped over the match information.
```

```
BLASTPS = 0
```

```
Close #BOFile
```

```
Kill BP.Path & "\" & BFileName
```

```
Kill BP.Path & "\" & BOFileName
```

```
Exit Function
```

```
BLASTPS = 0
```

```
E: Debug.Print "Error in BLASTPS"
```

```
Err.Raise Err.Number, Err.Description
```

```
End Function
```

```
Attribute VB_Name = "ColumnValidation"  
Option Explicit
```

```
Public Sub ColumnsValidate(RSPS As Recordset)
```

```
    'Function
```

```
    ' Check the validity of columns.
```

```
    'Arguments  
    ' RSPS: The parameter recordset for the Probesets being checked.
```

```
    'Notes
```

1. A column is valid if the values of the parameters used to calculate it are equal to the values of those parameters in the current instance of the appropriate parameter class.
2. This calculation affects settings in the global state variable pscolumns. It should therefore only be called on the Probesets and Probes used to build dbprobes, i.e. datasets.
3. Filter columns are always valid, since they are initialized with True.

```
Dim CreatePSValid As Boolean  
Dim LengthFilterValid As Boolean  
Dim PosFilterValid As Boolean  
Dim GCFilterValid As Boolean  
Dim dGDValid As Boolean  
Dim TMValid As Boolean  
Dim dGHValid As Boolean  
Dim dGMValid As Boolean  
Dim ClampValid As Boolean  
Dim RunValid As Boolean  
Dim HomologyValid As Boolean  
Dim Col As cColumn
```

```
'On startup, no columns exist, so bail.  
If pscolumns.Count = 0 Then Exit Sub
```

```
'If the recordsets are empty, then nothing is valid.
```

```
If Not IsGoodRS(RSPS) Then
```

```
    For Each Col In pscolumns
```

```
        If Col.CanBeInvalid Then Col.IsValid = False
```

```
    Next Col
```

```
    Exit Sub
```

```
End If
```

```
'set everything valid.
```

```
CreatePSValid = True
```

```
LengthFilterValid = True
```

```
PosFilterValid = True
```

```
GCFilterValid = True
```

```
dGDValid = True
```

```
TMValid = True
```

```
dGHValid = True
```

```
dGMValid = True
```

```
RunValid = True
```

```
ClampValid = True
```

```
HomologyValid = True
```

```
'Try to invalidate.
```

```
RSPS.MoveFirst
```

```
Do While Not RSPS.EOF
```

```
    CreatePSValid = CreatePSValid And CreatePSValid
```

```
    LengthFilterValid = LengthFilterValid And LengthFilterValid
```

```
    PosFilterValid = PosFilterValid And PosFilterValid
```

```
    GCFilterValid = GCFilterValid And GCFilterValid
```

```
    dGDValid = dGDValid And dGDValid
```

```
    TMValid = TMValid And TMValid
```

```
    dGHValid = dGHValid And dGHValid
```

```
    dGMValid = dGMValid And dGMValid
```

```
    RunValid = RunValid And RunValid
```

```
    ClampValid = ClampValid And ClampValid
```

```
    HomologyValid = HomologyValid And HomologyValid
```

```

dgmValid = dgmValid And dgmParams.Validate(RSPS) And
dgmFilterParams.Validate(RSPS)
RunValid = RunValid And RunParams.Validate(RSPS) And
RunFilterParams.Validate(RSPS)
ClampValid = ClampValid And ClampParams.Validate(RSPS) And
ClampFilterParams.Validate(RSPS)
HomologyValid = HomologyValid And HomologyParams.Validate(RSPS) And
HomologyFilterParams.Validate(RSPS)
RSPS.MoveNext
Loop

```

```

'Update columns.
PSColumns("Sequence").IsValid = CreatePSValid
PSColumns("Length").IsValid = CreatePSValid And LengthFilterValid
PSColumns("Length Filter").IsValid = PSColumns("Length").IsValid
PSColumns("Position").IsValid = CreatePSValid And PosFilterValid
PSColumns("Pos Filter").IsValid = PSColumns("Position").IsValid
PSColumns("GC").IsValid = CreatePSValid And GCFilterValid
PSColumns("GC Filter").IsValid = PSColumns("GC").IsValid
PSColumns("Duplex dG").IsValid = dgmValid
PSColumns("dG Filter").IsValid = PSColumns("Duplex dG").IsValid
PSColumns("TM").IsValid = TMValid
PSColumns("TM Filter").IsValid = PSColumns("TM").IsValid
PSColumns("dGH").IsValid = dgmValid
PSColumns("dGH Filter").IsValid = PSColumns("dGH").IsValid
PSColumns("dGM").IsValid = dgmValid
PSColumns("dGM Filter").IsValid = PSColumns("dGM").IsValid
PSColumns("Run Length").IsValid = RunValid
PSColumns("Run Filter").IsValid = PSColumns("Run Length").IsValid
PSColumns("Clamp").IsValid = ClampValid
PSColumns("Clamp Filter").IsValid = PSColumns("Clamp").IsValid
PSColumns("Homology").IsValid = HomologyValid
PSColumns("Homology Filter").IsValid = PSColumns("Homology").IsValid

```

'Set column heading backcolor according to validity.

```

For Each Col In PSColumns
    If Col.CanBeInValid Then
        If Col.IsValid Then
            frmMain.dbgProbes.Columns(Col.Name).ForeColor = vbBlack
        Else
            frmMain.dbgProbes.Columns(Col.Name).ForeColor = vbRed
        End If
    End If
Next Col

```

```

End Sub
Public Sub ColumnsExist(RSPS As Recordset, RSProbes As Recordset)
'-----

```

'Function
' Check the existence of columns.

'Arguments
' RSPS: The probeset(s) parameter recordset.
' RSProbes: The probeset probes.
'Notes

1. A column exists if
 - A. For each probeset, the appropriate existence parameter is set.

This indicates that some calculation has been performed to fill in the column, at sometime in the past. However, it is possible that

the selection criteria have changed since the calculation, bringing into scope probes that have not had the calculation performed. Thus, the existence parameter is necessary, but not sufficient.
B. For each probe, the appropriate value is not NULL. A lack of null values is necessary and sufficient; however, it is potentially time-consuming to check, so condition A is checked first, and condition B is checked only if A passes.

2. This calculation affects settings in the global state variable PSColumns. It should therefore only be called on the probesets and probes used to build dbgProbes, i.e. datSelps and datSelQuery.
3. Filter columns always exist, since they are initialized with True.

```

Dim Exists As Boolean
Dim Col As cColumn

```

```

'On startup, no columns exist, so bail.
If PSColumns.Count = 0 Then Exit Sub

```

'If the recordsets are empty, then nothing exists, nothing is valid.

```

If Not (IsGoodRS(RSPS) And IsGoodRS(RSProbes)) Then
    For Each Col In PSColumns
        If Col.CanUnExist Then Col.Exists = False
    Next Col
    Exit Sub
End If

```

'Probe creation calculations.

```

Exists = True
RSPS.MoveFirst
Do While Not RSPS.EOF
    Exists = Exists And RSPS("CreatePS-Exists")
    RSPS.MoveNext
Loop

```

```

PSColumns("Accession").Exists = Exists
PSColumns("PSName").Exists = Exists
PSColumns("Sequence").Exists = Exists
PSColumns("Length").Exists = Exists
PSColumns("Position").Exists = Exists
PSColumns("GC").Exists = Exists

```

'TM Calculations.

```

Exists = True
RSPS.MoveFirst
Do While Not RSPS.EOF
    Exists = Exists And RSPS("TM-Exists")
    RSPS.MoveNext
Loop

```

If Exists = True Then

```

    RSProbes.FindFirst "TM = NULL"
    If Not RSProbes.NoMatch Then Exists = False
End If
PSColumns("TM").Exists = Exists

```

'dGH Calculations.

```

Exists = True
RSPS.MoveFirst
Do While Not RSPS.EOF

```

```

Exists = Exists And RSPs("dgm-Exists")
RSPs.MoveNext
Loop
If Exists = True Then
    RSPs.FindFirst "dgm = NULL"
    If Not RSPs.NoMatch Then Exists = False
End If
PSColumns("dgm").Exists = Exists

'dgm Calculations.
Exists = True
RSPs.MoveFirst
Do While Not RSPs.EOF
    Exists = Exists And RSPs("dgm-Exists")
    RSPs.MoveNext
Loop
If Exists = True Then
    RSPs.FindFirst "dgm = NULL"
    If Not RSPs.NoMatch Then Exists = False
End If
PSColumns("dgm").Exists = Exists

'Clamp Calculations.
Exists = True
RSPs.MoveFirst
Do While Not RSPs.EOF
    Exists = Exists And RSPs("Clamp-Exists")
    RSPs.MoveNext
Loop
If Exists = True Then
    RSPs.FindFirst "[Clamp] = NULL"
    If Not RSPs.NoMatch Then Exists = False
End If
PSColumns("Clamp").Exists = Exists

'Run Length Calculations.
Exists = True
RSPs.MoveFirst
Do While Not RSPs.EOF
    Exists = Exists And RSPs("Run-Exists")
    RSPs.MoveNext
Loop
If Exists = True Then
    RSPs.FindFirst "[Run Length] = NULL"
    If Not RSPs.NoMatch Then Exists = False
End If
PSColumns("Run Length").Exists = Exists

'Homology Calculations.
Exists = True
RSPs.MoveFirst
Do While Not RSPs.EOF
    Exists = Exists And RSPs("Homology-Exists")
    RSPs.MoveNext
Loop
If Exists = True Then
    RSPs.FindFirst "[Homology] = NULL"
    If Not RSPs.NoMatch Then Exists = False

```

```

End If
PSColumns("Homology").Exists = Exists
End Sub

Attribute VB_Name = "CreateProbeset"
Option Explicit

Public Function CreatePSEqualL(Seq$, CPSP As cCreatePSPars, Probes$, -
    Positions$, GC#())
    -----
'Function:
' Create a set of probes using the EqualLength method.
'Parameters:
' Seq: The sequence from which the probes are created.
' CPSP: An instance of the parameter class for this algorithm.
' Probes: The array in which the probes will be returned.
' Positions: The array in which the starting positions will be returned.
' GC: The array in which GC content will be returned.
'Returns:
' The number of probes created.
'Notes:
' 2. The calling routine is responsible for adequately dimensioning Probes,
' Positions, and GC.
    -----
Dim StartPos% 'the starting position of a probe
Dim ProbeNum% 'the current probe

StartPos = 1
ProbeNum = LBound(Probes)
Do While StartPos <= Len(Seq) - CPSP.Length + 1
    If ProbeNum - Progress.StopAt = 0 Then Progress.CheckProgress StartPos
    Probes(ProbeNum) = Mid(Seq, StartPos, CPSP.Length)
    Positions(ProbeNum) = StartPos
    GC(ProbeNum) = DNA_GCContent(Probes(ProbeNum))
    StartPos = StartPos + CPSP.Spacing
    ProbeNum = ProbeNum + 1
Loop
CreatePSEqualL = ProbeNum - 1
End Function

Public Function CreatePSEqualTM(Seq$, CPSP As cCreatePSPars, Probes$, -
    Positions$, GC#())
    -----
'Function:
' Create a set of probes using the EqualTM method.
'Parameters:
' Seq: The sequence from which the probes are created.
' CPSP: An instance of the parameter class for this algorithm.
' Probes: The array in which the probes will be returned.
' Positions: The array in which the starting positions will be returned.
' GC: The array in which GC content will be returned.
'Returns:
' The number of probes created.
'Notes:
' 2. The calling routine is responsible for adequately dimensioning Probes,

```

```

' Positions, and GC.
'-----
Dim StartPos% 'the starting position of a probe
Dim ProbeNum% 'the current probe
Dim Length% 'current length of the probe being created
Dim TM# 'current TM of the probe being created

StartPos = 1
ProbeNum = LBound(Probes)
Do While StartPos <= Len(Seq)
    'Check progress
    If ProbeNum - Progress.Stopt = 0 Then Progress.CheckProgress StartPos
    'Create the probes, extending until melting point is above target.
    Length = 1
    Probes(ProbeNum) = Mid(Seq, StartPos, Length)
    Select Case CPSP.Duplex
        Case "DNA/DNA"
            TM = DNA_CalcTM(Probes(ProbeNum), CPSP.Conc)
        Case "DNA/RNA"
            TM = DR_CalcTM(Probes(ProbeNum), CPSP.Conc)
        End Select
    Length = Length + 1
    Do While TM < CPSP.Target
        If StartPos + Length > Len(Seq) Then
            ProbeNum = ProbeNum - 1
            Exit Do
        End If
        Probes(ProbeNum) = Mid(Seq, StartPos, Length)
        Select Case CPSP.Duplex
            Case "DNA/DNA"
                TM = DNA_CalcTM(Probes(ProbeNum), CPSP.Conc)
            Case "DNA/RNA"
                TM = DR_CalcTM(Probes(ProbeNum), CPSP.Conc)
        End Select
    Length = Length + 1
Loop
Positions(ProbeNum) = StartPos
GC(ProbeNum) = DNA_GCContent(Probes(ProbeNum))
ProbeNum = ProbeNum + 1
StartPos = StartPos + CPSP.Spacing
Loop
CreateSeqQualTM = ProbeNum - 1
End Function

Attribute VB_Name = "DNA_Manipulations"
Option Explicit

Public Function DNA_GCContent$(ByVal Seq$)
'Function
' Compute the GC content of a sequence.
'-----
Dim i%
Seq = LCase(Seq)
DNA_GCContent = 0
For i = 1 To Len(Seq)

```

```

If (Mid(Seq, i, 1) = "g" Or Mid(Seq, i, 1) = "c") Then DNA_GCContent =
DNA_GCContent + 1
Next i
DNA_GCContent = DNA_GCContent / Len(Seq) * 100
End Function

Public Function DNA_SeqTrim$(Seq$, Bases$, WhichEnd$)
'Function
' Trim off one end or the other of the sequence.
'-----
If Bases < Len(Seq) Then
    If WhichEnd = "5'" Then DNA_SeqTrim = Right(Seq, Len(Seq) - Bases)
    If WhichEnd = "3'" Then DNA_SeqTrim = Left(Seq, Len(Seq) - Bases)
End If
End Function

Public Function DNA_SeqKeep$(Seq$, Bases$, WhichEnd$)
'Function
' Remove all but the specified number of bases.
'-----
If Bases <= Len(Seq) Then
    If WhichEnd = "5'" Then DNA_SeqKeep = Left(Seq, Bases)
    If WhichEnd = "3'" Then DNA_SeqKeep = Right(Seq, Bases)
End If
End Function

Public Function DNA_RevComp$(ByVal Seq$)
'Function
' Given 5'-3' strand, returns 5'-3' complement.
'Revision
' 28-Jul-97: From PKWs routine of the same name.
'-----
Dim i%
Seq = LCase(Seq)
For i = Len(Seq) To 1 Step -1
    Select Case Mid$(Seq, i, 1)
        Case "a"
            DNA_RevComp = DNA_RevComp & "t"
        Case "c"
            DNA_RevComp = DNA_RevComp & "g"
        Case "g"
            DNA_RevComp = DNA_RevComp & "c"
        Case "t"
            DNA_RevComp = DNA_RevComp & "a"
        Case "u"
            DNA_RevComp = DNA_RevComp & "a"
        Case "n"
            DNA_RevComp = DNA_RevComp & "n"
        Case Else
            MsgBox "Unknown base found while calculating DNA_RevComp."
            DNA_RevComp = Seq
            Exit Function
    End Select
Next i
End Function

```

```
Public Function DNA_Comp$(ByVal Seq$)
```

```
Function
    Given 5'-3' strand, returns complement. This might be useless!
```

```
Dim i%
```

```
Seq = LCase(Seq)
```

```
For i = 1 To Len(Seq)
```

```
    Select Case Mid$(Seq, i, 1)
```

```
        Case "a"
```

```
            DNA_Comp = DNA_Comp & "t"
```

```
        Case "c"
```

```
            DNA_Comp = DNA_Comp & "g"
```

```
        Case "g"
```

```
            DNA_Comp = DNA_Comp & "c"
```

```
        Case "t"
```

```
            DNA_Comp = DNA_Comp & "a"
```

```
        Case "u"
```

```
            DNA_Comp = DNA_Comp & "a"
```

```
        Case "n"
```

```
            DNA_Comp = DNA_Comp & "n"
```

```
        Case Else
```

```
            MsgBox "Unknown base found while calculating DNA_Comp."
```

```
            DNA_Comp = Seq
```

```
        Exit Function
```

```
    End Select
```

```
Next i
```

```
End Function
```

```
Public Function DNA_Rev$(ByVal Seq$)
```

```
Function
```

```
    Given 5'-3' strand, returns 3'-5' sequence. This might be useless!
```

```
Dim i%
```

```
Seq = LCase(Seq)
```

```
For i = Len(Seq) To 1 Step -1
```

```
    Select Case Mid$(Seq, i, 1)
```

```
        Case "a"
```

```
            DNA_Rev = DNA_Rev & "a"
```

```
        Case "c"
```

```
            DNA_Rev = DNA_Rev & "c"
```

```
        Case "g"
```

```
            DNA_Rev = DNA_Rev & "g"
```

```
        Case "t"
```

```
            DNA_Rev = DNA_Rev & "t"
```

```
        Case "u"
```

```
            DNA_Rev = DNA_Rev & "u"
```

```
        Case "n"
```

```
            DNA_Rev = DNA_Rev & "n"
```

```
        Case Else
```

```
            MsgBox "Unknown base found while calculating DNA_Rev."
```

```
            DNA_Rev = Seq
```

```
        Exit Function
```

```
    End Select
```

```
Next i
```

```
End Function
```

```
Public Sub DNA_Str2Num(ByVal Seq$, ByRef NumSeq$(1))
```

```
Function
```

```
    Return a numeric array representing the DNA string.
```

```
Arguments
```

```
    Seq: The original sequence.
```

```
    NumSeq: An integer array representing the sequence, w/ lower bound 0.
```

```
Notes:
```

```
    1. NumSeq must be correctly dimensioned by the calling routine.
```

```
    since VB won't let me return an array from a function.
```

```
Dim i%
```

```
Seq = LCase(Seq)
```

```
For i = 0 To Len(Seq) - 1
```

```
    Select Case Mid$(Seq, i + 1, 1)
```

```
        Case "a"
```

```
            NumSeq(i) = 0
```

```
        Case "c"
```

```
            NumSeq(i) = 1
```

```
        Case "g"
```

```
            NumSeq(i) = 2
```

```
        Case "t"
```

```
            NumSeq(i) = 3
```

```
        Case "u"
```

```
            NumSeq(i) = 3
```

```
        Case "n"
```

```
            NumSeq(i) = -1
```

```
        Case Else
```

```
            MsgBox "Unknown base found while calculating DNA_Str2Num."
```

```
    End Select
```

```
Next i
```

```
End Sub
```

```
Public Function DNA_Num2Str$(NumSeq$(1))
```

```
Function:
```

```
    Return the string representation of the DNA string
```

```
Arguments:
```

```
    NumSeq: An integer array representing the sequence.
```

```
Dim i%
```

```
For i = 0 To UBound(NumSeq)
```

```
    Select Case NumSeq(i)
```

```
        Case "0"
```

```
            DNA_Num2Str = DNA_Num2Str & "a"
```

```
        Case "1"
```

```
            DNA_Num2Str = DNA_Num2Str & "c"
```

```
        Case "2"
```

```
            DNA_Num2Str = DNA_Num2Str & "g"
```

```
        Case "3"
```

```
            DNA_Num2Str = DNA_Num2Str & "t"
```

```
        Case "-1"
```

```
            DNA_Num2Str = DNA_Num2Str & "n"
```

```
        Case Else
```

```
            MsgBox "Unknown numeric code found while calculating DNA_Num2Str."
```

```
    End Select
```

```
Next i
```


End Function

Attribute VB Name = "Engine"

Option Explicit

Private Sub PSBlast(RSPS As Recordset, RSBLAST As Recordset)

Function

Perform all database gets/puts, etc, for BLAST calculation.

Arguments

RSPS: The parameters recordset (set parameters used in current record).

On Error GoTo E

Dim PSName\$

Dim RSBLAST As Recordset

Dim MatchAccessions\$(MaxBLASTHeaders)

Dim MatchHeaders\$(MaxBLASTHeaders)

Dim MatchScores\$(MaxBLASTHeaders)

Dim H\$

Dim MatchExpected\$(MaxBLASTHeaders)

Dim NumMatches\$

Dim FindStr\$

Get the probe set name and sequence.

PSName = RSPS("PSName")

Sequence = RSPS("CreatePS-Sequence")

If Len(Sequence) = 0 Then GoTo UpdateOnly

Run it.

Progress.ShowProgress "Calculating BLAST matches for ProbeSet = " & PSName, 0, 1

NumMatches = BLASTPS(PSName, Sequence, BlastPars, MatchAccession, MatchHeaders,

MatchScores, MatchExpected)

Update the results.

For H = 0 To NumMatches - 1

With RSBLAST

FindFirst "PSName = " & PSName & " AND Accession = " & MatchAccession(H) & "

If .NoMatch Then

AddNew

Edit

End If

Fields("PSName") = PSName

Fields("Accession") = MatchAccession(H)

Fields("Header") = MatchHeaders(H)

Fields("Score") = MatchScores(H)

Fields("Expect") = MatchExpected(H)

If MatchScores(H) > BlastPars.SCutoff Or MatchExpected(H) <

BlastPars.ECutoff Then

Fields("Use") = False

Else

Fields("Use") = True

End If

Update

End With

Next H

UpdateOnly:

Store parameters in the database.

BlastPars.StoreDB RSPS

Exit Sub

E: Debug.Print "Error in PSBlast"

Err.Raise Err.Number, , Err.Description

End Sub

Private Sub SeqCalcCreatePS(RSseq As Recordset, RSProbes As Recordset, -

RSPS As Recordset)

Function

Create a ProbeSet

Arguments

RSseq: The sequence recordset (positioned to current sequence).

RSPS: The Probes recordset (used to add a new ProbeSet).

RSProbes: The Probes recordset (used to add new probes).

RSPS: The Parameters recordset (used to add a new parameter record).

Errors Raised

ppErrorNoProbesCreated: Raised if no probes are generated.

Errors Handled

None.

Notes

1. Under normal use (e.g. when called by the GUI), this routine uses the

values in CreatePSpars.

2. Only RSseq must be a good recordset, since all the others are

used only to add records.

Dim Accession\$

Dim Seq\$

Dim SeqLength\$

Dim NumPS\$

Dim PSName\$

Dim Probes\$()

Dim Positions\$()

Dim GC\$()

Dim NumProbes\$

Dim P\$

Dim ProgressStr\$

On Error GoTo E

Check the recordset

If Not IsGoodRS(RSseq) Then Exit Sub

Get the sequence, accession, etc.

Accession = RSseq("Accession")

Seq = RSseq("Sequence")

SeqLength = RSseq("Length")

NumPS = RSseq("NumPS")

New probeset name.

```
PSName = Accession & " " & PSName & CStr(NumPS)
```

```
'Create temporary store for new probes.
Redim Probes (SeqLength)
Redim Positions (SeqLength)
Redim GC (SeqLength)
```

```
'Open the progress bar.
```

```
Progress.ShowProgress "Creating ProbeSet " & PSName, 0, NumProbes
```

```
'Create the probes.
```

```
Select Case CreatePSPar.Method
```

```
Case "Equal"
```

```
NumProbes = CreatePSEqual(Seq, CreatePSPar, Probes, Positions, GC)
```

```
Case "EqualTM"
```

```
NumProbes = CreatePSEqualTM(Seq, CreatePSPar, Probes, Positions, GC)
```

```
End Select
```

```
If NumProbes = 0 Then Err.Raise ppErrNoProbesCreated
```

```
Redim Preserve Probes (NumProbes)
```

```
Redim Preserve Positions (NumProbes)
```

```
Redim Preserve GC (NumProbes)
```

```
'Add a new ProbeSet.
```

```
RSSeq.Edit
```

```
RSSeq.Fields("NumPS") = NumPS + 1
```

```
RSSeq.Update
```

```
CreatePSPar.AddNewDB RSPar, Accession, PSName, SeqLength, Seq
```

```
'Create new records in the Probes table.
```

```
Progress.ShowProgress "Updating Probes in Database for ProbeSet " & PSName, 0,
```

```
NumProbes
```

```
For P = 0 To UBound(Probes) - 1
```

```
If P - Progress.StopAt = 0 Then Progress.CheckProgress P
```

```
RSProbes.AddNew
```

```
RSProbes.Fields("Accession") = Accession
```

```
RSProbes.Fields("PSName") = PSName
```

```
RSProbes.Fields("Sequence") = Probes(P)
```

```
RSProbes.Fields("Length") = Len(Probes(P))
```

```
RSProbes.Fields("Position") = Positions(P)
```

```
RSProbes.Fields("GC") = GC(P)
```

```
RSProbes.Update
```

```
Next P
```

```
'Clean exit.
```

```
Exit Sub
```

```
E: Debug.Print "Error in SeqCalcCreatePS"
```

```
Err.Raise Err.Number, Err.Description
```

```
End Sub
```

```
Public Sub SeqCalcEngine(Calcs)
```

```
-----
```

```
'Function
```

```
' SeqCalcEngine acts as the dispatcher for calculations
```

```
' performed on sequences. Each of the known calculations
```

```
' can be requested by specifying its name and subname.
```

```
' SeqCalcEngine then retrieves necessary records from the
```

```
' DB, then calls the specified calculation for each probe.
```

```
'Arguments
```

```
' Calc: The calculation to perform.
```

```
'Notes
```

```
' 1. The parameters used for a particular calculation are
```

```
' either made available in the optional parameter array
```

```
' Pairs, or retrieved from the class properties.
```

```
' 2. This routine loops over all selected sequences.
```

```
'Errors
```

```
' Errors are raised by not handled. Calling routine is responsible
```

```
-----
```

```
Dim RSSeqs As Recordset
```

```
Dim RSP As Recordset
```

```
Dim RSAllProbes As Recordset
```

```
Dim RSPS As Recordset
```

```
Dim PSName$
```

```
Dim ErrNumber
```

```
Dim ErrDescription
```

```
callee's
```

```
'Error handling.
```

```
On Error GoTo E
```

```
'Check that some sequences are available.
```

```
Set RSSeqs = frmMain.datSeqs.Recordset
```

```
If Not IsGoodRS(RSSeqs) Then Exit Sub
```

```
'Map recordset pointers.
```

```
Set RSP = frmMain.datP.Recordset
```

```
Set RSPS = frmMain.datPS.Recordset
```

```
'Close the grid.
```

```
frmMain.dbgProbes.Close
```

```
'Process all selected sequences.
```

```
RSSeqs.MoveFirst
```

```
Do While Not RSSeqs.EOF
```

```
'Begin a transaction block
```

```
'XXXXX -- causes an error when entering this routine for the second time, no
```

```
'idea why.
```

```
'BeginTrans
```

```
Select Case Calc
```

```
Case "CreatePS"
```

```
'Create the new probeSet.
```

```
SeqCalcCreatePS RSSeqs, RSP, RSPS
```

```
'Get PS Name
```

```
PSName = RSPS("PSName")
```

```
'Access all probes in this ProbeSet.
```

```
With frmMain.datOnePS
```

```
.RecordSource = "SELECT * FROM Probes WHERE PSName = '" & PSName &
```

```
""
```

```
.Refresh
```

```
Set RSAllProbes = .Recordset
```

```

End With

'Compute length and position filters.
If LengthFilterParms.Method <> "None" Then
    PSCalcLengthFilter RSAllProbes, RSFS
    PSColLengthFilter ("Length Filter").IsFilter = LengthFilterParms.AppImmediate
End If
If PosFilterParms.Method <> "None" Then
    PSCalcPosFilter RSAllProbes, RSFS
    PSColPosFilter ("Pos Filter").IsFilter = PosFilterParms.AppImmediate
End If
If GCFilterParms.Method <> "None" Then
    PSCalcGCFilter RSAllProbes, RSFS
    PSColGCFilter ("GC Filter").IsFilter = GCFilterParms.AppImmediate
End If
PSColLengthFilter ("Length").IsVisible = True
PSColPosFilter ("Position").IsVisible = True
PSColGCFilter ("GC").IsVisible = True

Case "3' Trim"
    SeqCalcTrim frmMain.datSelSeqs.Recordset, "3'"
Case "5' Trim"
    SeqCalcTrim frmMain.datSelSeqs.Recordset, "5'"
Case "3' Keep"
    SeqCalcKeep frmMain.datSelSeqs.Recordset, "3'"
Case "5' Keep"
    SeqCalcKeep frmMain.datSelSeqs.Recordset, "5'"

End Select

'End the transaction.
'CommitTrans

'Move to next selected sequence.
RS.SelSeqs.MoveNext

Loop

'Update appearances.
frmMain.Form_ChangeSequence
frmMain.Form_ChangePSName

'Reopen the grid.
frmMain.dbqProbes.ReOpen
Progress.Hide

'Clean exit.
Exit Sub

'Handle errors.
E: Debug.Print "Error in seqCalcEngine."
ErrNumber = Err.Number
ErrDescription = Err.Description

```

```

'Rollback
frmMain.Form_ChangeSequence
frmMain.Form_ChangePSName
frmMain.dbqProbes.ReOpen
Progress.Hide
frmMain.MousePointer = vbDefault
Err.Raise ErrNumber, , ErrDescription
End Sub

Private Sub SeqCalcTrim(RS As Recordset, WhichEnds)
'-----
'Function
' Trim bases from the 3' or 5' ends.
'-----
If Not IsGoodRS(RS) Then Exit Sub
RS.Edit
If WhichEnd = "5'" Then
    RS("Sequence") = DNA_SeqTrim(RS("sequence"), TrimParms.FivePTrim, "5'")
Else
    RS("Sequence") = DNA_SeqTrim(RS("sequence"), TrimParms.ThreePTrim, "3'")
End If
RS("Length") = Len(RS("sequence"))
RS("Accession") = RS("Accession") & ""
RS.Update
End Sub

Private Sub SeqCalcKeep(RS As Recordset, WhichEnds)
'-----
'Function
' Keep bases on the 3' or 5' ends.
'-----
If Not IsGoodRS(RS) Then Exit Sub
RS.Edit
If WhichEnd = "5'" Then
    RS("Sequence") = DNA_SeqKeep(RS("sequence"), TrimParms.FivePTrim, "5'")
Else
    RS("Sequence") = DNA_SeqKeep(RS("sequence"), TrimParms.ThreePTrim, "3'")
End If
RS("Length") = Len(RS("sequence"))
RS("Accession") = RS("Accession") & ""
RS.Update
End Sub

Private Sub SetField(RS As Recordset, Fields, Value)
'-----
'Function
' Set all entries for one field in the recordset.
'-----
RS.MoveFirst
Do While Not RS.EOF
    RS.Edit
    RS.Fields(Fields) = Value
    RS.Update
    RS.MoveNext
Loop
End Sub

```

```

Private Sub PutField(RS As Recordset, Fields, Value)
'-----
'Function
' Put all entries for one field in the recordset.
'Notes
' Since this operation can be timeconsuming, the progressbar is updated.
' So a progressbar has to be posted first!
'-----

```

```

Dim P%
P = 0
RS.MoveFirst
Do While Not RS.EOF
If P - Progress.StopAt = 0 Then Progress.CheckProgress P
RS.Edit
RS.Fields(Field) = Value(P)
RS.Update
P = P + 1
RS.MoveNext
Loop
End Sub

```

```

Private Sub GetField(RS As Recordset, Fields, Values)
'-----
'Function
' Get all entries for one field in the recordset.
'-----

```

```

Dim P%
P = 0
RS.MoveFirst
Do While Not RS.EOF
Values(P) = RS.Fields(Field)
RS.MoveNext
P = P + 1
Loop
End Sub

```

```

Private Sub PSCalcDGH(RSProbes As Recordset, RSPS As Recordset)
'-----

```

```

'Function
' Perform all database gets/puts, etc, for dGH calculation.
'Arguments
' RSProbes: The probes recordset (calculate dGH for all records).
' RSPS: The parameters recordset (set parameters used in current record).
'-----
On Error GoTo E

```

```

Dim PSName$
Dim NumProbes%
Dim Seq$()
Dim dGH#()

```

```

'Determine the number of probes.
NumProbes = NumRecords(RSProbes)
If NumProbes = 0 Then GoTo UpdateOnly

```

```

'Start it up.
PSName = RSProbes("PSName")
Progress.ShowProgress "Calculating dGH for Probeset " & PSName, 0, NumProbes

```

```

'Create space for database fields.
ReDim Seq(0 To NumProbes - 1)
ReDim dGH(0 To NumProbes - 1)

```

```

'Get the sequences.
GetField RSProbes, "Sequence", Seq

'Calculate dGH.
Select Case dGHPars.DR
Case "DNA"
DNA_CalcDGH Seq, dGHPars, dGH
Case "RNA"
RNA_CalcDGH Seq, dGHPars, dGH
End Select

```

```

'Update the results.
Progress.ShowProgress "Updating dGH in Database for Probeset " & PSName, 0,
NumProbes
PutField RSProbes, "dGH", dGH

```

```

UpdateOnly:

```

```

'Store parameters in the database, and update column visibility.
dGHPars.StoreDB RSPS

```

```

PSColumns("dGH").Exists = True
PSColumns("dGH").IsVisible = True

```

```

'Calculate filter
PSCalcDGHFilter RSProbes, RSPS

```

```

Exit Sub
E: Debug.Print "Error in PSCalcDGH"
Err.Raise Err.Number, , Err.Description
End Sub

```

```

Private Sub PSCalcDGM(RSProbes As Recordset, RSPS As Recordset)
'-----

```

```

'Function
' Perform all database gets/puts, etc, for dGM calculation.
'Arguments
' RSProbes: The probes recordset (calculate dGM for all records).
' RSPS: The parameters recordset (set parameters used in current record).
'-----

```

```

On Error GoTo E

```

```

Dim PSName$
Dim NumProbes%
Dim Seq$()
Dim dGM#()

```

```

'Determine the number of probes.
NumProbes = NumRecords(RSProbes)
If NumProbes = 0 Then GoTo UpdateOnly

```

```

'Start it up.
PSName = RSProbes("PSName")

```

```

Progress.ShowProgress "Calculating dGM for Probeset " & PSName, 0, NumProbes

'Create space for database fields.
ReDim Seq(0 To NumProbes - 1)
ReDim dGM(0 To NumProbes - 1)

'Get the sequences.
GetField RSProbes, "Sequence", Seq

'Calculate dGM.
Select Case dGMPars.DR
Case "DNA"
    DNA.CalcldGM Seq, dGMPars, dGM
Case "RNA"
    RNA.CalcldGM Seq, dGMPars, dGM
End Select

'Update the results.
Progress.ShowProgress "Updating dGM in Database for Probeset " & PSName, 0,
NumProbes
PutField RSProbes, "dGM", dGM

UpdateOnly:

'Store parameters in the database, and update column visibility.
dGMPars.StoreDB RSPS
PSColumns("dGM").Exists = True
PSColumns("dGM").IsVisible = True

'Calculate filter.
PSCalcldGMFilter RSProbes, RSPS

Exit Sub
E: Debug.Print "Error in PSCalcldGM."
Err.Raise Err.Number, , Err.Description
End Sub

```

```

Private Sub PSCalcRun(RSProbes As Recordset, RSPS As Recordset)
'Function
'Perform all database gets/puts, etc, for Run calculation.
'Arguments
'RSProbes: The probes recordset (calculate Run for all records).
'RSPS: The parameters recordset (set parameters used in current record).
-----
On Error GoTo E

Dim PSName$
Dim NumProbes$
Dim Pos$()
Dim Runs()

'Name of current probeset
'Number of probes in the recordset
'Position column from database.
'Run column, to be put to database.

'Determine the number of probes.
NumProbes = NumRecords(RSProbes)
If NumProbes = 0 Then GoTo UpdateOnly

'Update the results.
Progress.ShowProgress "Updating Run in Database for Probeset " & PSName, 0,
NumProbes
PutField RSProbes, "Run Length", Run

```

```

'Start it up.
PSName = RSProbes("PSName")
Progress.ShowProgress "Calculating Run for Probeset " & PSName, 0, NumProbes

'Create space for database fields.
ReDim Pos(0 To NumProbes - 1)
ReDim Run(0 To NumProbes - 1)

'Get the sequences.
GetField RSProbes, "Position", Pos

'Calculate Run.
CalcRun Pos, RunPars, Run

'Update the results.
Progress.ShowProgress "Updating Run in Database for Probeset " & PSName, 0,
NumProbes
PutField RSProbes, "Run Length", Run

UpdateOnly:

'Store parameters in the database, and update column visibility.
RunPars.StoreDB RSPS
PSColumns("Run Length").Exists = True
PSColumns("Run Length").IsVisible = True

'Calculate the run filter.
PSCalcRunFilter RSProbes, RSPS

Exit Sub
E: Debug.Print "Error in PSCalcRun"
Err.Raise Err.Number, , Err.Description
End Sub

Private Sub PSCalcHomology(RSProbes As Recordset, RSPS As Recordset, RSBLast As
Recordset)
'Function
'Perform all database gets/puts, etc, for Homology calculation.
'Arguments
'RSProbes: The probes recordset (calculate Homology for all records).
'RSPS: The parameters recordset (set parameters used in current record).
-----
On Error GoTo E

Dim PSName$
Dim NumProbes$
Dim Seq$()
Dim Homology$()

'Name of current probeset
'Number of probes in the recordset
'Sequences of probes
'Run column, to be put to database

'Determine the number of probes.
NumProbes = NumRecords(RSProbes)
If NumProbes = 0 Then GoTo UpdateOnly

'Start it up.
'Note that the progress bar in this case will show progress over sequences
'in the database. Since we don't know the total, we'll set the maximum to

```

```

'100, and have the called routine do fraction calculations.
PSName = RSProbes("PSName")
Progress.ShowProgress "Calculating Homology for Probeset " & PSName, 0, 100

'Create space for database fields.
Redim Seq(0 To NumProbes - 1)
Redim Homology(0 To NumProbes - 1)

'Get the sequences.
GetField RSProbes, "Sequence", Seq

'Calculate Homology.
CalcHomology Seq, HomologyParms, Homology, RSblast, PSName

'Update the results.
Progress.ShowProgress "Updating Homology in Database for Probeset " & PSName, 0,
NumProbes
PutField RSProbes, "Homology", Homology

UpdateOnly:

'Store parameters in the database, and update column visibility.
HomologyParms.StoreDB RSps
HomologyParms("Homology").Exists = True
PSColumns("Homology").IsVisible = True

'calculate the homology filter.
PSCalcHomofilter RSProbes, RSps

Exit Sub
E: Debug.Print "Error in PSCalcHomology"
Err.Raise Err.Number, , Err.Description
End Sub

Private Sub PSCalcTM(RSProbes As Recordset, RSps As Recordset)
'-----
'Function
' Perform all database gets/puts, etc, for TM calculation.
'Arguments
' RSProbes: The probes recordset (calculate TM for all records).
' RSps: The parameters recordset (set parameters used in current record).
'-----
On Error GoTo E

Dim PSName$ 'name of current probeset
Dim NumProbes$ 'number of probes in the recordset
Dim Seq$() 'sequence column from database.
Dim TM$() 'TM column, to be put to database.

'Determine the number of probes.
NumProbes = NumRecords(RSProbes)
If NumProbes = 0 Then GoTo UpdateOnly

'Start it up.
PSName = RSProbes("PSName")
Progress.ShowProgress "Calculating TM for Probeset " & PSName, 0, NumProbes

'Create space for database fields.

```

```

Redim Seq(0 To NumProbes - 1)
Redim TM(0 To NumProbes - 1)

'Get the sequences.
GetField RSProbes, "Sequence", Seq

'Calculate TM.
Select Case TMPars.Duplex
Case "DNA/DNA"
DNA.CalcAllTM Seq, TMPars, TM
Case "DNA/RNA"
DR.CalcAllTM Seq, TMPars, TM
End Select

'Update the results.
Progress.ShowProgress "Updating TM in Database for Probeset " & PSName, 0,
NumProbes
PutField RSProbes, "TM", TM

UpdateOnly:

'Store parameters in the database, and update column visibility.
TMPars.StoreDB RSps
PSColumns("TM").Exists = True
PSColumns("TM").IsVisible = True

'Calculate the filter.
PSCalcTMfilter RSProbes, RSps

Exit Sub
E: Debug.Print "Error in PSCalcTM"
Err.Raise Err.Number, , Err.Description
End Sub

Private Sub PSCalcDGD(RSProbes As Recordset, RSps As Recordset)
'-----
'Function
' Perform all database gets/puts, etc, for dGD calculation.
'Arguments
' RSProbes: The probes recordset (calculate dGD for all records).
' RSps: The parameters recordset (set parameters used in current record).
'-----
On Error GoTo E

Dim PSName$ 'name of current probeset
Dim NumProbes$ 'number of probes in the recordset
Dim Seq$() 'sequence column from database.
Dim dGD$() 'dGD column, to be put to database.

'Determine the number of probes.
NumProbes = NumRecords(RSProbes)
If NumProbes = 0 Then GoTo UpdateOnly

'Start it up.
PSName = RSProbes("PSName")
Progress.ShowProgress "Calculating dG (Duplex) for Probeset " & PSName, 0,
NumProbes

```

```
'Create space for database fields.
Redim Seq(0 To NumProbes - 1)
Redim dGd(0 To NumProbes - 1)
```

```
'Get the sequences.
GetField RSProbes, "Sequence", Seq
```

```
'Calculate dGd.
Select Case dGdPars.Duplex
Case "DNA"
DNA_CalcldGd Seq, dGdPars, dGd
Case "RNA"
DR_CalcldGd Seq, dGdPars, dGd
End Select
```

```
'Update the results.
Progress.ShowProgress "Updating dG (Duplex) in Database for Probeset " & PSName,
0, NumProbes
PutField RSProbes, "Duplex dG", dGd
```

```
UpdateOnly:
```

```
'Store parameters in the database, and update column visibility.
dGdPars.StoreDB RSFS
PSColumns("Duplex dG").Exists = True
PSColumns("Duplex dG").IsVisible = True
```

```
'Calculate the filter.
```

```
PSCalcldGdFilter RSProbes, RSFS
```

```
Exit Sub
```

```
E: Debug.Print "Error in PSCalcldGd"
```

```
Err.Raise Err.Number, , Err.Description
```

```
End Sub
```

```
Private Sub PSCalcClamp(RSProbes As Recordset, RSFS As Recordset)
```

```
Function
```

```
Perform all database gets/puts, etc, for Clamp calculation.
```

```
Arguments
```

```
RSProbes: The probes recordset (calculate Clamp for all records).
```

```
RSFS: The parameters recordset (set parameters used in current record).
```

```
On Error GoTo E
```

```
Dim PSNames
```

```
'name of current probeset
```

```
Dim NumProbes
```

```
'number of probes in the recordset
```

```
Dim Seqs()
```

```
'sequence column from database.
```

```
Dim Clamp#()
```

```
'Clamp column, to be put to database.
```

```
'Determine the number of probes.
```

```
NumProbes = NumRecords(RSProbes)
```

```
If NumProbes = 0 Then GoTo UpdateOnly
```

```
'Start it up.
```

```
PSName = RSProbes("PSName")
```

```
Progress.ShowProgress "Calculating Clamp for Probeset " & PSName, 0, NumProbes
```

```
© Copyright Hewlett-Packard Company, 1998 29 Attorney Docket No. 10971464-1
```

```
'Create space for database fields.
Redim Seq(0 To NumProbes - 1)
Redim Clamp(0 To NumProbes - 1)
```

```
'Get the sequences.
```

```
GetField RSProbes, "Sequence", Seq
```

```
'Calculate Clamp.
```

```
Select Case ClampPars.Duplex
```

```
Case "DNA/DNA"
```

```
DNA_CalcldClamp Seq, ClampPars, Clamp
```

```
Case "DNA/RNA"
```

```
'XXXXX
```

```
DNA_CalcldClamp Seq, ClampPars, Clamp
```

```
End Select
```

```
'Update the results.
```

```
Progress.ShowProgress "Updating Clamp in Database for Probeset " & PSName, 0,
NumProbes
```

```
PutField RSProbes, "Clamp", Clamp
```

```
UpdateOnly:
```

```
'Store parameters in the database, and update column visibility.
```

```
ClampPars.StoreDB RSFS
```

```
PSColumns("Clamp").Exists = True
```

```
PSColumns("Clamp").IsVisible = True
```

```
'Calculate the filter
```

```
PSCalcldClampFilter RSProbes, RSFS
```

```
Exit Sub
```

```
E: Debug.Print "Error in PSCalcldClamp"
```

```
Err.Raise Err.Number, , Err.Description
```

```
End Sub
```

```
Public Sub PSCalcEngine(Calc#)
```

```
Function
```

```
PSCalcEngine acts as the dispatcher for calculations
```

```
performed on Probesets. Each of the known calculations
```

```
can be requested by specifying its name.
```

```
Arguments
```

```
Calc: The calculation to perform.
```

```
Notes
```

```
1. All numeric calculations proceed by
```

```
- checking that the sequences column exists.
```

```
- checking that the current parameters set for this calculation
```

```
are the same as those previously used, and, if not, nulling
```

```
out all previous calculations.
```

```
- performing the calculation.
```

```
2. Filter calculations proceed as above, with the preliminary step
```

```
of checking the existence of the column to be filtered. If it isn't
```

```
there, the calculation is done.
```

```
3. The Probes grid is closed before DB updates, then opened
```

```
when they are all over, to reduce screen thrashing.
```

```
© Copyright Hewlett-Packard Company, 1998 30 Attorney Docket No. 10971464-1
```



```

If Not dGHFilterParms.Validate(RSPS) Then SetField RSAllProbes, "dGH
Filter", True
PSCalcCDGH RSQProbes, RSPS

Case "dGH Filter"
If Not PSColumns("dGH").Exists Then
If dGHPars.Exists(RSPS) And Not dGHPars.Validate(RSPS) Then SetField
RSAllProbes, "dGH", True
If Not dGHFilterParms.Validate(RSPS) Then SetField RSAllProbes, "dGH
Filter", True
Else
PSCalcCDGH RSQProbes, RSPS
If Not dGHFilterParms.Validate(RSPS) Then SetField RSAllProbes, "dGH
Filter", True
PSCalcCDGHFilter RSQProbes, RSPS
End If

Case "dGM"
If dGMPars.Exists(RSPS) And Not dGMPars.Validate(RSPS) Then SetField
RSAllProbes, "dGM", Null
If Not dGMFilterParms.Validate(RSPS) Then SetField RSAllProbes, "dGM
Filter", True
PSCalcCDGM RSQProbes, RSPS

Case "dGM Filter"
If Not PSColumns("dGM").Exists Then
If dGMPars.Exists(RSPS) And Not dGMPars.Validate(RSPS) Then SetField
RSAllProbes, "dGM", True
If Not dGMFilterParms.Validate(RSPS) Then SetField RSAllProbes, "dGM
Filter", True
Else
PSCalcCDGM RSQProbes, RSPS
If Not dGMFilterParms.Validate(RSPS) Then SetField RSAllProbes, "dGM
Filter", True
PSCalcCDGMFilter RSQProbes, RSPS
End If

Case "Run"
If RunPars.Exists(RSPS) And Not RunPars.Validate(RSPS) Then SetField
RSAllProbes, "Run Length", Null
If Not RunFilterParms.Validate(RSPS) Then SetField RSAllProbes, "Run
Filter", True
PSCalcRun RSQProbes, RSPS

Case "Run Filter"
If Not PSColumns("Run Length").Exists Then
If RunPars.Exists(RSPS) And Not RunPars.Validate(RSPS) Then SetField
RSAllProbes, "Run", True
If Not RunFilterParms.Validate(RSPS) Then SetField RSAllProbes, "Run
Filter", True
PSCalcRun RSQProbes, RSPS
Else
If Not RunFilterParms.Validate(RSPS) Then SetField RSAllProbes, "Run
Filter", True
PSCalcRunFilter RSQProbes, RSPS
End If

Case "Clamp"

```

```

If ClampPars.Exists(RSPS) And Not ClampPars.Validate(RSPS) Then SetField RSAllProbes, "Clamp", Null
If Not ClampFilterPars.Validate(RSPS) Then SetField RSAllProbes, "Clamp Filter", True
PSCalcClamp RSQProbes, RSPS

Case "Clamp Filter"
If Not PSColumns("Clamp").Exists Then
If ClampPars.Exists(RSPS) And Not ClampPars.Validate(RSPS) Then SetField RSAllProbes, "Clamp", Null
If Not ClampFilterPars.Validate(RSPS) Then SetField RSAllProbes, "Clamp Filter", True
PSCalcClamp RSQProbes, RSPS
Else
If Not ClampFilterPars.Validate(RSPS) Then SetField RSAllProbes, "Clamp Filter", True
PSCalcClampFilter RSQProbes, RSPS
End If

Case "Length Filter"
If Not PSColumns("Length").Exists Then Exit Sub
If Not LengthFilterPars.Validate(RSPS) Then SetField RSAllProbes, "Length Filter", True
PSCalcLengthFilter RSQProbes, RSPS

Case "Pos Filter"
If Not PSColumns("Position").Exists Then Exit Sub
If Not PosFilterPars.Validate(RSPS) Then SetField RSAllProbes, "Pos Filter", True
PSCalcPosFilter RSQProbes, RSPS

Case "GC Filter"
If Not PSColumns("GC").Exists Then Exit Sub
If Not GCFilterPars.Validate(RSPS) Then SetField RSAllProbes, "GC Filter", True
PSCalcGCFilter RSQProbes, RSPS

Case "BLAST"
PSCalcBlast RSPS, RSblast

Case "Homology"
If HomologyPars.Exists(RSPS) And Not HomologyPars.Validate(RSPS) Then SetField RSAllProbes, "Homology", Null
If Not HomofilterPars.Validate(RSPS) Then SetField RSAllProbes, "Homology Filter", True
PSCalcHomology RSQProbes, RSPS, RSblast

Case "Homology Filter"
If Not PSColumns("Homology").Exists Then
If HomologyPars.Exists(RSPS) And Not HomologyPars.Validate(RSPS) Then SetField RSAllProbes, "Homology", Null
If Not HomofilterPars.Validate(RSPS) Then SetField RSAllProbes, "Homology Filter", True
PSCalcHomology RSQProbes, RSPS, RSblast
Else
If Not HomofilterPars.Validate(RSPS) Then SetField RSAllProbes, "Homology Filter", True
PSCalcHomofilter RSQProbes, RSPS

```

```

End If
End Select

'End the transaction.
CommitTrans

'Move to next selected probeset.
RSselp.NextMoveNext

Loop

'Reopen the grid.
Progress.Hide
frmMain.dbgProbes.ReOpen
DoEvents

'Update appearances.
frmMain.frm_ChangePSName

'Clean exit.
frmMain.MousePointer = vbDefault
Exit Sub

'Handle errors.
E: Debug.Print "Error in PSCalcEngine."
ErrNumber = Err.Number
ErrDescription = Err.Description
Rollback
frmMain.frm_ChangePSName
frmMain.dbgProbes.ReOpen
Progress.Hide
frmMain.MousePointer = vbDefault
Err.Raise ErrNumber, ErrDescription

End Sub

```

```

Private Sub PSCalcTMFilter(RSPROBES As Recordset, RSPS As Recordset)
'Function
' Perform all database gets/puts, etc, for TM filter calculation.
'Arguments
' RSPROBES: The probes recordset (calculate TM Filter for all records).
' RSPS: The parameter recordset (set parameters used in current record).
-----
On Error GoTo E

Dim PSName$
Dim NumProbes
Dim TM#()
Dim TMFilter() As Boolean
'Name of current probeset
'number of probes in the recordset
'TM column from database
'TM Filter column to database

'Determine the number of probes.
NumProbes = NumRecords(RSPROBES)
If NumProbes = 0 Then GoTo UpdateOnly

'Create space for database fields.
ReDim dgdTM(0 To NumProbes - 1)

```

```

'Start it up.
PSName = RSPROBES("PSName")
Progress.ShowProgress "Calculating TM Filter for Probeset " & PSName, 0,
NumProbes

'Create space for database fields.
ReDim TM(0 To NumProbes - 1)
ReDim TMFilter(0 To NumProbes - 1)

'Get the TMs.
GetField RSPROBES, "TM", TM

'Calculate TM Filter.
CalcTMFilter TM, TMFilterParams, TMFilter

'Update the results.
Progress.ShowProgress "Updating TM Filter in Database for Probeset " & PSName,
0, NumProbes
PutField RSPROBES, "TM Filter", TMFilter

UpdateOnly:
'Store parameters in the database, and update column visibility.
TMFilterParams.StoreDB RSPS
PSCOLUMNS("TM Filter").IsFilter = TMFilterParams.ApplyImmediate

Exit Sub
E: Debug.Print "Error in PSCalcTMFilter"
Err.Raise Err.Number, Err.Description
End Sub

Private Sub PSCalcDGDFilter(RSPROBES As Recordset, RSPS As Recordset)
'Function
' Perform all database gets/puts, etc, for dGD filter calculation.
'Arguments
' RSPROBES: The probes recordset (calculate dGD Filter for all records).
' RSPS: The parameter recordset (set parameters used in current record).
-----
On Error GoTo E

Dim PSName$
Dim NumProbes
Dim dgd#()
Dim dgdFilter() As Boolean
'Name of current probeset
'number of probes in the recordset
'dGD column from database
'dGD Filter column to database

'Determine the number of probes.
NumProbes = NumRecords(RSPROBES)
If NumProbes = 0 Then GoTo UpdateOnly

'Start it up.
PSName = RSPROBES("PSName")
Progress.ShowProgress "Calculating dG (Duplex) Filter for Probeset " & PSName,
0, NumProbes

'Create space for database fields.
ReDim dgd(0 To NumProbes - 1)
ReDim dgdFilter(0 To NumProbes - 1)

```

```

'Get the dGDs.
GetField RSProbes, "Duplex dG", dGD

'Calculate dGD Filter.
CalcGDFilter dGD, dGDFilterPars, dGDFilter

'Update the results.
Progress.ShowProgress "Updating dG (Duplex) Filter in Database for Probeset " &
PSName, 0, NumProbes
PutField RSProbes, "dGD Filter", dGDFilter

UpdateOnly:

'store parameters in the database, and update column visibility.
dGDFilterPars.StoreDB RSPS
PSColumns("dGD Filter").IsFilter = dGDFilterPars.AppImmediate

Exit Sub
E: Debug.Print "Error in PSCalcGDFilter"
Err.Raise Err.Number, , Err.Description
End Sub

Private Sub PSCalcClampFilter(RSPRObes As Recordset, RSPS As Recordset)
'-----
'Function
' Perform all database gets/puts, etc, for Clamp filter calculation.
' Arguments
' RSPRObes: The probes recordset (calculate Clamp Filter for all records).
' RSPS: The parameter recordset (set parameters used in current record).
'-----
On Error Goto E

Dim PSName$
Dim NumProbes$
Dim Clamp#()
Dim ClampFilter() As Boolean
Dim ClampFilterColumn$

'Name of current probeset
'number of probes in the recordset
'dGH column from database
'dGH Filter column to database

'Determine the number of probes.
NumProbes = NumRecords(RSPRObes)
If NumProbes = 0 Then Goto UpdateOnly

'Start it up.
PSName = RSPRObes("PSName")
Progress.ShowProgress "Calculating Clamp Filter for Probeset " & PSName, 0,
NumProbes

'Create space for database fields.
ReDim Clamp(0 To NumProbes - 1)
ReDim ClampFilter(0 To NumProbes - 1)

'Get the dGDs.
GetField RSProbes, "dGD", dGD

'Calculate dGD Filter.
CalcGDFilter dGD, dGDFilterPars, dGDFilter

'Update the results.
Progress.ShowProgress "Updating dG Filter in Database for Probeset " & PSName,
0, NumProbes
PutField RSProbes, "dGD Filter", dGDFilter

UpdateOnly:

'store parameters in the database, and update column visibility.
dGDFilterPars.StoreDB RSPS

```

```

'Update the results.
Progress.ShowProgress "Updating Clamp Filter in Database for Probeset " &
PSName, 0, NumProbes
PutField RSProbes, "Clamp Filter", ClampFilter

UpdateOnly:

'store parameters in the database, and update column visibility.
ClampFilterPars.StoreDB RSPS
PSColumns("Clamp Filter").IsFilter = ClampFilterPars.AppImmediate

Exit Sub
E: Debug.Print "Error in PSCalcClampFilter"
Err.Raise Err.Number, , Err.Description
End Sub

Private Sub PSCalcdGHFilter(RSPRObes As Recordset, RSPS As Recordset)
'-----
'Function
' Perform all database gets/puts, etc, for dGH filter calculation.
' Arguments
' RSPRObes: The probes recordset (calculate dGH Filter for all records).
' RSPS: The parameter recordset (set parameters used in current record).
'-----
On Error Goto E

Dim PSName$
Dim NumProbes$
Dim dGH#()
Dim dGHFilter() As Boolean
Dim dGHFilterColumn$

'Name of current probeset
'number of probes in the recordset
'dGH column from database
'dGH Filter column to database

'Determine the number of probes.
NumProbes = NumRecords(RSPRObes)
If NumProbes = 0 Then Goto UpdateOnly

'Start it up.
PSName = RSPRObes("PSName")
Progress.ShowProgress "Calculating dGH Filter for Probeset " & PSName, 0,
NumProbes

'Create space for database fields.
ReDim dGH(0 To NumProbes - 1)
ReDim dGHFilter(0 To NumProbes - 1)

'Get the dGHs.
GetField RSProbes, "dGH", dGH

'Calculate dGH Filter.
CalcdGHFilter dGH, dGHFilterPars, dGHFilter

'Update the results.
Progress.ShowProgress "Updating dGH Filter in Database for Probeset " & PSName,
0, NumProbes
PutField RSProbes, "dGH Filter", dGHFilter

UpdateOnly:

'store parameters in the database, and update column visibility.
dGHFilterPars.StoreDB RSPS

```

PSCOLUMNS("dgm Filter").IsFilter = dgmFilterPars.AppImmediate

Exit Sub
E: Debug.Print "Error in PSCalcDGMFilter"
Err.Raise Err.Number, , Err.Description
End Sub

Private Sub PSCalcGCFilter(RSPROBES As Recordset, RSPS As Recordset)

Function
Perform all database gets/puts, etc, for GC filter calculation.
Arguments
ASPROBES: The probes recordset (calculate GC Filter for all records).
RSPS: The parameter recordset (set parameters used in current record).
On Error GOTO E

Dim PSName\$ 'name of current probeset
Dim NumProbes\$ 'number of probes in the recordset
Dim GC#() 'GC column from database
Dim GCFilter() As Boolean 'GC Filter column to database

Determine the number of probes.
NumProbes = NumRecords(RSPROBES)
If NumProbes = 0 Then Goto UpdateOnly

Start it up.
PSName = RSPROBES("PSName")
Progress.ShowProgress "Calculating GC Filter for Probeset " & PSName, 0,
NumProbes

Create space for database fields.
ReDim GC(0 To NumProbes - 1)
ReDim GCFilter(0 To NumProbes - 1)

Get the GCs.
GetField RSPROBES, "GC", GC

Calculate GC Filter.
CalcGCFilter GC, GCFilterPars, GCFilter

Update the results.
Progress.ShowProgress "Updating GC Filter in Database for Probeset " & PSName,
0, NumProbes
PutField RSPROBES, "GC Filter", GCFilter

UpdateOnly:

Store parameters in the database, and update column visibility.
GCFilterPars.StoredB RSPS
PSCOLUMNS("GC Filter").IsFilter = GCFilterPars.AppImmediate

Exit Sub
E: Debug.Print "Error in PSCalcGCFilter"
Err.Raise Err.Number, , Err.Description
End Sub

Private Sub PSCalcDGMFilter(RSPROBES As Recordset, RSPS As Recordset)

Function
Perform all database gets/puts, etc, for dgm filter calculation.
Arguments
RSPROBES: The probes recordset (calculate dgm Filter for all records).
RSPS: The parameter recordset (set parameters used in current record).
On Error GOTO E

Dim PSName\$ 'name of current probeset
Dim NumProbes\$ 'number of probes in the recordset
Dim dgm#() 'dgm column from database
Dim dgmFilter() As Boolean 'dgm Filter column to database

Determine the number of probes.
NumProbes = NumRecords(RSPROBES)
If NumProbes = 0 Then Goto UpdateOnly

Start it up.
PSName = RSPROBES("PSName")
Progress.ShowProgress "Calculating dgm Filter for Probeset " & PSName, 0,
NumProbes

Create space for database fields.
ReDim dgm(0 To NumProbes - 1)
ReDim dgmFilter(0 To NumProbes - 1)

Get the dgm's.
GetField RSPROBES, "dgm", dgm

Calculate dgm Filter.
CalcDGMFilter dgm, dgmFilterPars, dgmFilter

Update the results.
Progress.ShowProgress "Updating dgm Filter in Database for Probeset " & PSName,
0, NumProbes
PutField RSPROBES, "dgm Filter", dgmFilter

UpdateOnly:

Store parameters in the database, and update column visibility.
dgmFilterPars.StoredB RSPS
PSCOLUMNS("dgm Filter").IsFilter = dgmFilterPars.AppImmediate

Exit Sub
E: Debug.Print "Error in PSCalcDGMFilter"
Err.Raise Err.Number, , Err.Description
End Sub

Private Sub PSCalcRunFilter(RSPROBES As Recordset, RSPS As Recordset)

Function
Perform all database gets/puts, etc, for Run filter calculation.
Arguments
RSPROBES: The probes recordset (calculate Run Filter for all records).
RSPS: The parameter recordset (set parameters used in current record).
On Error GOTO E

```

Dim PSName$
'Name of current probeset
'Dim NumProbes%
'number of probes in the recordset
'Dim Run%()
'Run column from database
'Dim RunFilter() As Boolean
'Run Filter column to database

'Determine the number of probes.
NumProbes = NumRecords(RSProbes)
If NumProbes = 0 Then GoTo updateOnly

'Start it up.
PSName = RSProbes("PSName")
PProgress.ShowProgress "Calculating Run Filter for ProbeSet " & PSName, 0,
NumProbes

'Create space for database fields.
ReDim Run(0 To NumProbes - 1)
ReDim RunFilter(0 To NumProbes - 1)

'Get the Runs.
GetField RSProbes, "Run Length", Run

'Calculate Run Filter.
calcRunFilter Run, RunFilterPars, RunFilter

'Update the results.
PProgress.ShowProgress "Updating Run Filter in Database for ProbeSet " & PSName, 0, NumProbes
PutField RSProbes, "Run Filter", RunFilter

updateOnly:

'Store parameters in the database, and update column visibility.
RunFilterPars.storeDB RSPS
PSColumns("Run Filter").IsFilter = RunFilterPars.AppImmediate

Exit Sub
E: Debug.Print "Error in PSCalcRunFilter"
Err.Raise Err.Number, , Err.Description
End Sub

Private Sub PSCalcHomofilter(RSProbes As Recordset, RSPS As Recordset)
-----
'Function
' Perform all database gets/puts, etc, for Homology filter calculation.
'Arguments
' RSProbes: The probes recordset (calculate Homology Filter for all records).
' RSPS: The parameter recordset (set parameters used in current record).
-----
On Error GoTo E

Dim PSName$
'Name of current probeset
Dim NumProbes%
'number of probes in the recordset
Dim Homologys()
'Homology column from database
Dim Homofilter() As Boolean
'Homology Filter column to database

'Determine the number of probes.
NumProbes = NumRecords(RSProbes)

```

```

If NumProbes = 0 Then GoTo UpdateOnly

'start it up.
PSName = RSProbes("PSName")
Progress.ShowProgress "Calculating Homology Filter for ProbeSet " & PSName, 0, NumProbes

'Create space for database fields.
ReDim Homology(0 To NumProbes - 1)
ReDim HomoFilter(0 To NumProbes - 1)

'Get the Homologys.
GetField RSProbes, "Homology", Homology

'Calculate Homology Filter.
CalcHomoFilter Homology, HomoFilterPars, HomoFilter

'Update the results.
Progress.ShowProgress "Updating Homology Filter in Database for ProbeSet " & PSName, 0, NumProbes
PutField RSProbes, "Homology Filter", HomoFilter

UpdateOnly:

'Store parameters in the database, and update column visibility.
HomoFilterPars.StoreDB RSProbes
PFCColumns("Homology Filter").IsFilter = HomoFilterPars.AppImmediate

Exit Sub
E: Debug.Print "Error in PSCalcHomoFilter"
Err.Raise Err.Number, , Err.Description
End Sub

Private Sub PSCalcLengthFilter(RSProbes As Recordset, RSAs As Recordset)
'Function
'Perform all database gets/puts, etc, for Length filter calculation.
'Arguments
'RSProbes: The probes recordset (calculate Length Filter for all records).
'RSAs: The parameter recordset (set parameters used in current record).
-----
On Error GoTo E

Dim PSNames$
Dim NumProbes$
Dim Length$()
Dim LengthFilter() As Boolean

'Name of current probeset
'Number of probes in the recordset
'Length column from database
'Length Filter column to database

'Determine the number of probes.
NumProbes = NumRecords(RSProbes)
If NumProbes = 0 Then GoTo UpdateOnly

'start it up.
PSName = RSProbes("PSName")
Progress.ShowProgress "Calculating Length Filter for ProbeSet " & PSName, 0, NumProbes

```

```
'Create space for database fields.
Redim Length(0 To NumProbes - 1)
Redim LengthFilter(0 To NumProbes - 1)
```

```
'Get the Lengths.
GetField RSProbes, "Length", Length
```

```
'Calculate Length Filter.
CalcLengthFilter Length, LengthFilterPars, LengthFilter
```

```
'Update the results.
Progress.ShowProgress "Updating Length Filter in Database for ProbeSet " &
PSName, 0, NumProbes
PutField RSProbes, "Length Filter", LengthFilter
UpdateOnly:
```

```
'Store parameters in the database, and update column visibility.
LengthFilterPars.StoreDB RSPs
PSColumns("Length Filter").IsFilter = LengthFilterPars.AppImmediate
```

```
Exit Sub
E: Debug.Print "Error in PSCalcLengthFilter"
Err.Raise Err.Number, , Err.Description
End Sub
```

```
Private Sub PSCalcPosFilter(RSProbes As Recordset, RSPs As Recordset)
```

```
Function
    Perform all database gets/puts, etc, for Pos filter calculation.
```

```
Arguments
    RSProbes: The probes recordset (calculate Pos Filter for all records).
    RSPs: The parameter recordset (set parameters used in current record).
```

```
On Error GoTo E
```

```
Dim PSName$
Dim NumProbes$
Dim NumPos()
Dim PosFilter() As Boolean
Dim SeqLengths
    'name of current probeset
    'number of probes in the recordset
    'pos column from database
    'Pos Filter column to database
    'Length of the originating sequence.
```

```
'Determine the number of probes.
NumProbes = NumRecords(RSProbes)
If NumProbes = 0 Then Goto UpdateOnly
```

```
'Start it up.
PSName = RSProbes("PSName")
Progress.ShowProgress "Calculating Pos Filter for ProbeSet " & PSName, 0,
NumProbes
```

```
'Create space for database fields.
Redim Pos(0 To NumProbes - 1)
Redim PosFilter(0 To NumProbes - 1)
```

```
'Get the Positions.
GetField RSProbes, "Position", Pos
```

```
'Get the sequence length.
SeqLength = RSPs("CreatePS-SeqLength")
```

```
'Calculate Pos Filter.
CalcPosFilter Pos, SeqLength, PosFilterPars, PosFilter
```

```
'Update the results.
Progress.ShowProgress "Updating Pos Filter in Database for ProbeSet " & PSName,
0, NumProbes
PutField RSProbes, "Pos Filter", PosFilter
UpdateOnly:
```

```
'Store parameters in the database, and update column visibility.
PosFilterPars.StoreDB RSPs
PSColumns("Pos Filter").IsFilter = PosFilterPars.AppImmediate
```

```
Exit Sub
E: Debug.Print "Error in PSCalcPosFilter"
Err.Raise Err.Number, , Err.Description
End Sub
```

```
Attribute VB_Name = "Filters"
Option Explicit
```

```
Public Sub CalcTMFilter(TM() As TMFPars As CTMFilterPars, Filter() As Boolean)
```

```
Function
    Decide which probes pass the TM filter.
```

```
Arguments
    TM: An array of TMs
    TMFPars: An instance of the parameter class CTMFilterPars.
    Filter: The return array of filter output values.
```

```
'The TM filter can be implemented directly by the generic filter.
CalcGenericFilter TM, TMFPars, Filter
End Sub
```

```
Public Sub CalcGDGFilter(dGD() As dGDFPars As cdGDFilterPars, Filter() As
Boolean)
```

```
Function
    Decide which probes pass the dGD filter.
```

```
Arguments
    dGD: An array of dGDs
    dGDFPars: An instance of the parameter class cdGDFilterPars.
    Filter: The return array of filter output values.
```

```
'The dGD filter can be implemented directly by the generic filter.
CalcGenericFilter dGD, dGDFPars, Filter
End Sub
```

Public Sub CalcGenericFilter(Values, FilterPars As Object, Filter() As Boolean)

```

'Function
'  Decide which values pass the filter.
'Arguments
'  Values: An array of values to be filtered against.
'  FilterPars: The filter parameters
'  Filter: The returned filter settings.
'Notes
'  This routine is used to implement the common filter methods. FilterPars is
'  declared as an object rather than a specific class, and thus needs to have
'  only the members accessed on the execution path.
'  FilterPars must have data member Method.
'  If method is None, no other members are required.
'  If method is Min, then member Min is required.
'  If method is Max, then member Max is required.
'  If method is InRange, then members Max and Min are required.
'  If method is Distance, then members Target and Distance are required.
'  If method is Nbest, then members Target and Nbest are required.
'  If method is Percent, then members Target and Percent are required.
'History
'  31-Jul-97: PW

```

```

Dim Distance#()
Dim Target#
Dim MinValue#, MaxValue#
Dim Index#()
Dim NumValues#
Dim T#
Dim N#

'absolute difference between values and target
'target value
'min and max of value
'index array to track sort
'number of values we are working with
'index
'number to set.

```

'Determine the number of probes we are calculating filter for.

NumValues = UBound(Values) + 1

'If method is "None", set all to true.

If FilterPars.Method = "None" Then

For T = 0 To NumValues - 1

Filter(T) = True

Next T

Exit Sub

End If

'If method is "Min", set all >= min to true.

If FilterPars.Method = "Min" Then

For T = 0 To NumValues - 1

Filter(T) = (Values(T) >= Val(FilterPars.Min))

Next T

Exit Sub

End If

'If method is "Max", set all <= max to true.

If FilterPars.Method = "Max" Then

For T = 0 To NumValues - 1

Filter(T) = (Values(T) <= Val(FilterPars.Max))

Next T

Exit Sub

End If

```

'If method is "InRange", set all >= min and <= max to true.
If FilterPars.Method = "InRange" Then
For T = 0 To NumValues - 1
Filter(T) = ((Values(T) >= Val(FilterPars.Min)) And (Values(T) <=
Val(FilterPars.Max)))
Next T
Exit Sub
End If

```

'For other methods, we have to do more work!

'Size the arrays.

ReDim Distance(NumValues - 1)

ReDim Index(NumValues - 1)

'Find the minimum and maximum values.

For T = 0 To NumValues - 1

If MinValue > Values(T) Then MinValue = Values(T)

If MaxValue < Values(T) Then MaxValue = Values(T)

Next T

'Choose the target value.

If (FilterPars.Method = "Distance") Or _

(FilterPars.Method = "Nbest") Or _

(FilterPars.Method = "Percent") _

Then Target = FilterPars.Target

If (FilterPars.Method = "Nlowest") Or

(FilterPars.Method = "PercentLowest") _

Then Target = MinValue

If (FilterPars.Method = "NHighest") Or

(FilterPars.Method = "PercentHighest") _

Then Target = MaxValue

'Compute distances of values from target; set up index and filter.

For T = 0 To NumValues - 1

Distance(T) = Abs(Values(T) - Target)

Index(T) = T

Filter(T) = False

Next T

'If method is "Distance", set all elements within distance as true.

If FilterPars.Method = "Distance" Then

For T = 0 To NumValues - 1

Filter(T) = Distance(T) <= FilterPars.Distance

Next T

Exit Sub

End If

'Choose the number to set

If (FilterPars.Method = "Nbest") Or _

(FilterPars.Method = "Nlowest") Or _

(FilterPars.Method = "NHighest") _

Then N = FilterPars.Nbest

If (FilterPars.Method = "Percent") Or

(FilterPars.Method = "PercentLowest") Or _

(FilterPars.Method = "PercentHighest") _

Then N = NumValues * FilterPars.Percent / 100 - 1

```

If N >= NumValues Then N = NumValues - 1
'Sort on distance from target.
Quicksort Distance, Index, 0, NumValues - 1
'Set the best.
For T = 0 To N - 1
    Filter(Index(T)) = True
Next T
End Sub

```

```

Public Sub CalcRunFilter(Run#(), RunFPars As CRunFilterPars, Filter() As
Boolean)

```

```

'Function
' Decide which probes pass the Run filter.
'Arguments
' Run: An array of Runs
' RunFPars: An instance of the parameter class CRunFilterPars.
' Filter: The return array of filter output values.
Dim OldMins

```

```

'The special method "All" can be handled by using the method Min, with min=1.
If RunFPars.Method = "All" Then
    OldMin = RunFPars.Min
    RunFPars.Min = 1
    RunFPars.Method = "Min"
    CalcGenericFilter Run, RunFPars, Filter
    RunFPars.Method = "All"
    RunFPars.Min = OldMin
Exit Sub
End If

```

```

'Otherwise, filters are generic.
CalcGenericFilter Run, RunFPars, Filter
End Sub

```

```

Public Sub CalcdgmFilter(dgm#(), dgmFPars As cdgmFilterPars, Filter() As
Boolean)

```

```

'Function
' Decide which probes pass the dgm filter.
'Arguments
' dgm: An array of dgms
' dgmFPars: An instance of the parameter class cdgmFilterPars.
' Filter: The return array of filter output values.

```

```

CalcGenericFilter dgm, dgmFPars, Filter
End Sub

```

```

Public Sub CalcClampFilter(Clamp#(), ClampFPars As cClampFilterPars, Filter() As
Boolean)

```

```

'Function
' Decide which probes pass the Clamp Filter.
'Arguments
' Clamp: An array of clamps
' ClampFPars: An instance of the parameter class cClampFilterPars.
' Filter: The return array of filter output values.
'Notes
'History
' 31-Jul-97: PW

```

```

CalcGenericFilter Clamp, ClampFPars, Filter
End Sub

```

```

Public Sub CalcGcFilter(GC#(), GCFPars As cGcFilterPars, Filter() As Boolean)

```

```

'Function
' Decide which probes pass the GC filter.
'Arguments
' GC: An array of GCs
' GCFPars: An instance of the parameter class cGcFilterPars.
' Filter: The return array of filter output values.

```

```

CalcGenericFilter GC, GCFPars, Filter
End Sub

```

```

Public Sub CalcdghFilter(dgh#(), dghFPars As cdghFilterPars, Filter() As
Boolean)

```

```

'Function
' Decide which probes pass the dgh filter.
'Arguments
' dgh: An array of dghs
' dghFPars: An instance of the parameter class cdghFilterPars.
' Filter: The return array of filter output values.

```

```

CalcGenericFilter dgh, dghFPars, Filter
End Sub

```

```

Public Sub CalcHomoFilter(Homology#(), HFPars As cHomoFilterPars, Filter() As
Boolean)

```

```

'Function
' Decide which probes pass the homology filter.
'Arguments
' Homology: An array of homologies
' HFPars: An instance of the parameter class cHomoFilterPars.
' Filter: The return array of filter output values.

```

```

CalcGenericFilter Homology, HFPars, Filter

```


End Sub

Public Sub CalcPosFilter(Pos() As Integer, SeqLength As Integer, PPars As cPosFilterPars, Filter() As Boolean)

Function
Decide which probes pass the position filter.

Arguments

Pos: An array of positions.

SeqLength: The length of the sequence from which the probeset was derived.

PPars: An instance of the parameter class cPosFilterPars.

Filter: The return array of filter output values.

The target for this filter should be the specified end.

If PPars.WhichEnd = "5-prime" Then

PPars.Target = 1

Else

PPars.Target = SeqLength

End If

CalcGenericFilter Pos, PPars, Filter

End Sub

Public Sub CalcLengthFilter(Length() As Integer, LFPars As cLengthFilterPars, Filter() As Boolean)

Function

Decide which probes pass the length filter.

Arguments

Length: An array of lengths.

LFPars: An instance of the parameter class cLengthFilterPars.

Filter: The return array of filter output values.

CalcGenericFilter length, LFPars, Filter

End Sub

Attribute VB_Name = "GenBank"

Option Explicit

Private Const GBTerminator\$ = "/"

Private Const MaxGBFileLength\$ = 100000 "maximum allowed GenBank file length"

Public Sub EntrezLoadDB(Accession\$, RS As Recordset, KeepFile As Boolean)

Function

Grab a sequence using the Entrez CGI.

Arguments

Accession: The accession of the required sequence.

RS: The recordset into which this sequence should be loaded.

KeepFile: Controls whether the downloaded file is retained or removed.

On Error GoTo E

Dim GBFileName\$

Dim EntrezData\$, EntrezData2\$ "GB Entries.

Dim i%

'Grab the data.
EntrezData = frmEntrezGrab.Entrez.GetGenBankData (Accession, "")

'Create CRIFs out of LFs.

For i = 1 To Len(EntrezData)

If Mid\$(EntrezData, i, 1) = vbLf Then EntrezData2 = EntrezData2 + vbCr

EntrezData2 = EntrezData2 + Mid\$(EntrezData, i, 1)

Next i

'Write the acquired data out to a file.

GBFileName = frmMain.cdgSeqs.InitDir & "\" & Accession & ".gb"

Open GBFileName For Output As #1

Print #1, EntrezData2

Close #1

'Load the database.

SeqLoadDB ReadGBRecord(GBFileName), RS

'Remove the file.

If Not KeepFile Then Kill GBFileName

Exit Sub

E: Debug.Print "Error in EntrezLoadDB"

Err.Raise Err.Number, Err.Description

End Sub

Public Sub SeqLoadDB(SeqRC As SeqRecord, RS As Recordset)

Function

Create a new record in the Sequence table, populate from SeqRC

Arguments

SeqRC: A sequence structure, as read from file, etc.

RS: The recordset into which the sequence should be loaded.

Errors:

Attempting to create a new sequence record whose Accession matches a

previously inserted sequence causes a DB error.

On Error GoTo E

BeginTrans

AllowDbgClick = False

RS.AddNew

RS.Fields("Header") = SeqRC.Header

RS.Fields("Accession") = SeqRC.Accession

RS.Fields("Locus") = SeqRC.Locus

RS.Fields("Length") = SeqRC.Length

RS.Fields("Sequence") = SeqRC.Sequence

RS.Fields("Selected") = True

RS.Update

RS.Bookmark = RS.LastModified

CommitTrans

DoEvents

AllowDbgClick = True

Exit Sub

E: Debug.Print "Error in SeqLoadDB"

Rollback

DoEvents

```

AllowdbgClick = True
Err.Raise Err.Number, , Err.Description
End Sub

Public Function ReadGBRecord(FileName) As SeqRecord
-----
'Function
' Read a GenBank record from a file into a SeqRecord.
' Arguments
'   FileName: The filename of the GB record.
'Returns
'   The GenBank record as a SeqRecord.
'Errors:
'   Returned to caller.
'   1. Check file is less than MaxGBFileLength
'   2. Check that all required fields are present.
-----
On Error GoTo E

Dim GBFile$
Dim FileLength$
Dim LineLength$
Dim GBText$
Dim GBItem$
Dim KeyField$
Dim foo$

Dim GB As SeqRecord
Dim the returned record

'open the file
GBFile = FreeFile
Open FileName For Binary As #GBFile

'read in the file as one long string
FileLength = FileLen(FileName)
If FileLength > MaxGBFileLength Then Err.Raise ppErrGBFileLength
GBText = Input(FileLength, #GBFile)

'check for required fields
If Instr(GBText, "LOCUS") = 0 Then Err.Raise ppErrGBFileFormat, , "No LOCUS
field found in GB file."
If Instr(GBText, "ACCESSION") = 0 Then Err.Raise ppErrGBFileFormat, , "No
ACCESSION field found in GB file."
If Instr(GBText, GBTerminator) = 0 Then Err.Raise ppErrGBFileFormat, , "No
terminator (//) found in GB file."

'process items
Do While True

'extract an item
LineLength = Instr(GBText, vbCrLf) - 1
GBItem = Left(GBText, LineLength)
GBHeader = GBHeader & GBItem & vbCrLf
GBText = Right(GBText, FileLength - LineLength - 2)
FileLength = FileLength - LineLength - 2

If Left(GBItem, 2) = GBTerminator Then

```

```

GoTo breakwhile
End If

'extract keyfield, and process item
KeyField = StrField(GBItem, LineLength, 10)
'if keyfield is not numeric -> header
If IsNumeric(KeyField) = False Then
'if keyfield occupies first column -> keyword
If Left(KeyField, 1) <> " " Then
'process keywords
Select Case KeyField
Case "LOCUS"
foo = StrField(GBItem, LineLength, 2)
GB.Locus = StrField(GBItem, LineLength, 10)
GB.Length = StrField(GBItem, LineLength, 7)
Case "ACCESSION"
foo = StrField(GBItem, LineLength, 2)
GB.Accession = StrField(GBItem, LineLength, 6)
End Select
'if keyfield has only two blanks -> subkeyword
Elseif Left(KeyField, 3) <> " " Then
'if keyfield has only 4 blanks -> feature code
Elseif Left(KeyField, 5) <> " " Then
'it must be a continuation line
Else
End If
Else
'feature has started
Do While LineLength > 10
GB.Sequence = GB.Sequence & Left(GBItem, 10)
GBItem = Right(GBItem, LineLength - 11)
LineLength = LineLength - 11
Loop
If LineLength > 0 Then
GB.Sequence = GB.Sequence & GBItem
End If
Loop
breakwhile:
Close GBFile

'Strip any whitespace.
GB.Sequence = PackSequence(GB.Sequence)

'check sequence length.
If GB.Length <> Len(GB.Sequence) Then
MsgBox "Sequence length differs from header specification."
GB.Length = Len(GB.Sequence)
GB.Accession = GB.Accession & ""
End If

ReadGBRecord = GB
Exit Function
E: Debug.Print "Error in ReadGBRecord"
Err.Raise Err.Number, , Err.Description
End Function

```

```
Public Function ReadFASTARecord(FASTAFile#) As SeqRecord
```

```

Function
    ' Read a FASTA record from a file into a SeqRecord.
    ' Arguments
    '   FASTA: The file descriptor for the FASTA file.
    ' Returns
    '   The FASTA record as a SeqRecord.
    ' Errors:
    '   Returned to caller.
    ' Notes:
    '   The FASTA file must already be opened for input, and should be positioned
    '   at the start of a FASTA header. It will read an entire sequence, leaving
    '   the file positioned on the subsequent header.
    '-----
    On Error GoTo E

    Dim LineTexts
    Dim SR As SeqRecord

    'one line of the file
    'the returned record

    'Read the header line, check it is legal.
    Line Input #FASTAFile, lineText
    If Left(lineText, 1) <> ">" Then Err.Raise ppErrFASTAFormat, "Header line does
    not begin w/ >."

```

```
"Extract the length.
LineText = Right(LineText, Len(LineText) - Instr(LineText, "=") - 3)
If Instr(LineText, "/") > 0 Then LineText = Left(LineText, Instr(LineText, "/"
- 1)
SR.Length = CLng(LineText)
```

```
'Read sequence lines till we are done.
Do While (Len(SR.Sequence) <> SR.Length)
  Line Input #FASTAfile, Linetext
  SR.Sequence = SR.Sequence & Linetext
Loop
```

ReadFASTRecord = SR

```
Exit Function
E: Debug.Print "Error in ReadFASTRecord"
Err.Raise Err.Number, Err.Description
End Function
```

```
public Function ReadFASTAHeader(FASTAFile f) As SeqRecord
```

```

Function
    Read a FASTA header from a file into a SeqRecord.

```

Arguments

FASTA: The file descriptor for the FASTA header file.

'Returns

• The FASTA record as

Errors

```

: Returned to caller.
:
: Notes:
: The FASTA header file must already be opened for input, and should be
: positioned
: at the start of a FASTA header.

```

On Error Goto E

```

Dim LineText$
'one line of the file
'line processed to find length
Dim LenText$
'line processed to find accession
Dim AccText$
'the returned record
Dim SR As SeqRecord

'Read the header line, check it is legal.
Line Input #FASTAFile, LineText
If Left(LineText, 1) <> ">" Then Err.Raise ppErrFASTAFormat, "Header line does
not begin w/ '>'"

'Extract the length.
LenText = Right(LineText, Len(LineText) - Instr(LineText, "len=") - 3)
If Instr(LineText, "#") <> 0 Then LenText = Left(LineText, Instr(LenText, "#") -
1)
SR.Length = CInt(LenText)

'Extract the accession.
AccText = Right(LineText, Len(LineText) - Instr(LineText, "gb=") - 2)
If Instr(AccText, "#") <> 0 Then AccText = Left(AccText, Instr(AccText, "#") -
1)
SR.Accession = LTrim(RTrim(AccText))

ReadFASTAHeader = SR

Exit Function
E: Debug.Print "Error in ReadFASTARecord"
Err.Raise Err.Number, Err.Description
End Function

```

Public Function GBCookSeq(ByVal RawSeq) As String
 'Convert raw sequence (no spaces, etc) into more
 'palatable form, with spaces every 10 nucleotides
 'and newlines every 60 nucleotides.

```

Dim Bases$ 'count the bases
Dim LineLength$ 'remaining lineLength

LineLength = Len(RawSeq)
Bases = 1
Do While LineLength > 60
    GBCookSeq = GBCookSeq & _
        Format(Format(Bases, "#####"), "#####") & " "
    Bases = Bases + 60
    Dim If
    For i = 1 To 6
        GBCookSeq = GBCookSeq & StrField(RawSeq, LineLength, 10) & " "
    Next i
    GBCookSeq = GBCookSeq & vbCrLf
Loop
If Len(RawSeq) <> 0 Then
    GBCookSeq = GBCookSeq & _
        Format(Format(Bases, "#####"), "#####") & " "
    Do While Len(RawSeq) > 10
        GBCookSeq = GBCookSeq & StrField(RawSeq, LineLength, 10) & " "
    
```

```

Loop
  GBCookSeq = GBCookSeq & RawSeq
End If
End Function

Public Function RTFGBCookSeq(yVal RawSeq) As String
  'Convert raw sequence (no spaces, etc) into more
  'palatable form, with spaces every 10 nucleotides
  'and newlines every 60 nucleotides.

  Dim Bases$
  'count the bases
  Dim LineLength$
  'remaining line length

  LineLength = Len(RawSeq)
  Bases = 1
  Do While LineLength > 60
    RTFGBCookSeq = RTFGBCookSeq &
    Format(Format(Bases, "#####"), "#####") & " "
    Bases = Bases + 60
    Dim i$
    For i = 1 To 6
      RTFGBCookSeq = RTFGBCookSeq & StrField(RawSeq, LineLength, 10) & " "
    Next i
    RTFGBCookSeq = RTFGBCookSeq & "\par"
  Loop
  If Len(RawSeq) <> 0 Then
    RTFGBCookSeq = RTFGBCookSeq &
    Format(Format(Bases, "#####"), "#####") & " "
    Do While Len(RawSeq) > 10
      RTFGBCookSeq = RTFGBCookSeq & StrField(RawSeq, LineLength, 10) & " "
    Loop
  End If
  RTFGBCookSeq = RTFGBCookSeq & RawSeq
End Function

Attribute VB_Name = "SQL"
Option Explicit

Public Function BuildPonePSQuery(PSName)
  'Function
  ' Build an SQL Query to select Probes from the Probes table
  ' that pass all filters, and come from the named ProbeSet.
  'Arguments
  ' Accession: The accession to search on.
  ' PSName: The ProbeSet name to search on.
  'Notes
  ' 1. Accession or PSName = <none> is allowed, and results
  ' in an empty recordset.
  ' 2. All fields corresponding to a column (that is, all fields
  ' named in the Columns table) are returned, to match the layout
  ' of the Probes grid. The visibility of various columns is
  ' controlled by the column Visible properties.
  ' 3. The currently set filters are used as a further selection
  ' on the probes returned.
  Dim SelectC$, FromC$, FilterC$, WhereC$, OrderByC$
  Dim Cs, Fc

```

```

'Set up the select clause -- return all columns.
SelectC = "SELECT *"
'For C = 1 To PSColumns.Count
' If C <> 1 Then SelectC = SelectC & ", "
' SelectC = SelectC & "[" & PSColumns(C).Name & "]"
'Next C

'Set up the the from clause -- always the Probes table.
FromC = "FROM Probes "

'Set up the where clause
For F = 1 To PSColumns.Count
  If (PSColumns(F).IsFilter = True) Then
    If FilterC = "" Then
      FilterC = "[" & PSColumns(F).Name & "]" = True) "
    Else
      FilterC = FilterC & " AND (" & PSColumns(F).Name & "]" = True) "
    End If
  End If
Next F
WhereC = "WHERE ( PSName = " & PSName & " )"
If PSColumns("User Filter").IsFilter Then
  If FilterC <> "" Then
    WhereC = WhereC & " AND ( (User Filter) = 2) OR ( (User Filter) = 1)
    AND " & FilterC & " ) "
  Else
    WhereC = WhereC & " AND ( (User Filter) <> 0) ) "
  End If
Else
  If FilterC <> "" Then
    WhereC = WhereC & " AND ( " & FilterC & " ) "
  Else
    WhereC = WhereC & " "
  End If
End If

'put it all together.
BuildPonePSQuery = SelectC & FromC & WhereC & " ";

End Function

Public Function BuildProbesQuery()
  'Function
  ' Build an SQL Query to select Probes from the Probes table.
  'Arguments
  ' Accession: The accession to search on.
  ' PSName: The ProbeSet name to search on.
  'Notes
  ' 1. Accession or PSName = <none> is allowed, and results
  ' in an empty recordset.
  ' 2. All fields corresponding to a column (that is, all fields
  ' named in the Columns table) are returned, to match the layout
  ' of the Probes grid. The visibility of various columns is
  ' controlled by the column Visible properties.
  ' 3. The currently set filters are used as a further selection
  ' on the probes returned.

```

```

'-----
Dim SelectC$, FromC$, FilterC$, WhereC$, OrderByC$
Dim C$, P$
Dim Col As cColumn

'Set up the select clause -- return all columns.
SelectC = "SELECT *"
For C = 1 To PSColumns.Count
    If C <> 1 Then SelectC = SelectC & ", "
    SelectC = SelectC & "Probes.[ " & PSColumns(C).Name & " ] "
Next C

'Set up the the from clause -- from the SelectProbes query.
FromC = "FROM SelectProbes "

'Set up the where clause -- filter the rows.
For F = 1 To PSColumns.Count
    For Each Col In PSColumns
        If (Col.IsFilter = True) Then
            If FilterC = "" Then
                FilterC = "[ " & Col.Name & " ] = True) "
            Else
                FilterC = FilterC & " AND ([ " & Col.Name & " ] = True) "
            End If
        End If
    Next Col
WhereC = "WHERE ( (Selected = True) "
If PSColumns("User Filter").IsFilter Then
    If FilterC <> "" Then
        WhereC = WhereC & " AND ( ([User Filter] = 1) "
        AND " & FilterC & " ) ) "
    Else
        WhereC = WhereC & " AND ([User Filter] <> 0) ) "
    End If
Else
    If FilterC <> "" Then
        WhereC = WhereC & " AND ( " & FilterC & " ) "
    Else
        WhereC = WhereC & " "
    End If
End If

'Set up the orderby clause.
'Default to sorting by position, then check columns.
OrderByC = "Position "
For Each Col In PSColumns
    If (Col.IsSort = "Ascending") Then OrderByC = "[ " & Col.Name & " ] ASC "
    If (Col.IsSort = "Descending") Then OrderByC = "[ " & Col.Name & " ] DESC "
Next Col
OrderByC = "ORDER BY Probes.Accession, Probes.PSName, " & OrderByC

'put it all together.
BuildProbesQuery = SelectC & FromC & WhereC & OrderByC & ";"

End Function

Attribute VB_Name = "Stats"

```

```

Option Explicit
Public Type StatRecord
    Count As Long
    Min As Double
    Ten As Double
    Median As Double
    Ninety As Double
    Max As Double
    Range As Double
    Average As Double
    Std As Double
End Type

Public Function Statistics(Vector#()) As StatRecord
'-----
'Function
' Calculate all statistics of a vector.
'-----
Dim L#, U#
Dim Min#, Max#, Total#, Total2#
Dim Index#()
Dim i#

'Get array bounds.
L = LBound(Vector)
U = UBound(Vector)

'Set the count.
Statistics.Count = UBound(Vector) - LBound(Vector) + 1
ReDim Index(L To U)

'Take a pass through to get min, max, total, total^2, index.
Min = Vector(L)
Max = Vector(U)
For i = L To U
    Index(i) = i
    Total = Total + Vector(i)
    Total2 = Total2 + Vector(i) * Vector(i)
    If Vector(i) > Max Then Max = Vector(i)
    If Vector(i) < Min Then Min = Vector(i)
Next i

'Set min, max, range, average, std.
Statistics.Min = Min
Statistics.Max = Max
Statistics.Range = Max - Min
Statistics.Average = Total / Statistics.Count
Statistics.Std = Sqr((Total2 / Statistics.Count) - Statistics.Average ^ 2)

'Sort.
QuickSort Vector, Index, L, U

'Percentiles.
Statistics.Ten = Vector(L + (Statistics.Count - 1) * 0.1)
Statistics.Median = Vector(L + (Statistics.Count - 1) * 0.5)
Statistics.Ninety = Vector(L + (Statistics.Count - 1) * 0.9)

```

End Function

Attribute VB_Name = "Thermo"
Option Explicit

'There is no good reason to have two sets of thermodynamic parameters,
'other than history...

'Parameters used by TM calculations.

Private Const Rgas# = 1.987
Private Const DNA_DuplexH#(0 To 3, 0 To 3)

enthalpies.

Private DNA_DuplexH#(0 To 3, 0 To 3)

Private DNA_InitGCS#, DNA_InitATS#

Private DNA_SelfIS#

Private DNA_EndTAH#

Private DR_DuplexH#(0 To 3, 0 To 3)

enthalpies.

Private DR_DuplexH#(0 To 3, 0 To 3)

Private DR_InitH#, DR_InitS#

enthalpy/entropy.

'Parameters used by hairpin calculations.

'Zuker provides H and G @ 37. So on load, calculate S, then recalculate G as
needed.

Private Const NumZukerLoops# = 30

Private ZLoopLengths#(0 To NumZukerLoops - 1)

Private Const DNA_ZmiscLoop# = 1.079

Private Const RNA_ZmiscLoop# = 1.079

entropies.

Private Const ZminLoops# = 4

Private Const NumTetraLoops# = 8

Private ZTetraLoops#(0 To NumTetraLoops - 1)

'DNA free energy.

Private DNAAHP_StackG#(0 To 3, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private DNAAHP_TStackH#(0 To 3, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private DNAAHP_DangleG#(0 To 1, 0 To 3, 0 To 3, 0 To 3)

Private DNAAHP_LoopG#(0 To 2, 0 To NumZukerLoops - 1)

Private DNAAHP_TloopG#(0 To NumTetraLoops - 1)

'DNA enthalpy.

Private DNAAHP_StackH#(0 To 3, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private DNAAHP_TStackH#(0 To 3, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private DNAAHP_DangleH#(0 To 1, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private DNAAHP_LoopH#(0 To 2, 0 To NumZukerLoops - 1)

Private DNAAHP_TloopH#(0 To NumTetraLoops - 1)

'DNA entropy.

Private DNAAHP_Stacks#(0 To 3, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private DNAAHP_Tstacks#(0 To 3, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private DNAAHP_Dangles#(0 To 1, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private DNAAHP_Loops#(0 To 2, 0 To NumZukerLoops - 1)

Private DNAAHP_Tloops#(0 To NumTetraLoops - 1)

'RNA free energy.

Private RNAHP_StackG#(0 To 3, 0 To 3, 0 To 3, 0 To 3, 0 To 3)
Private RNAHP_TStackH#(0 To 3, 0 To 3, 0 To 3, 0 To 3, 0 To 3)
Private RNAHP_DangleG#(0 To 1, 0 To 3, 0 To 3, 0 To 3, 0 To 3)
Private RNAHP_LoopG#(0 To 2, 0 To NumZukerLoops - 1)
Private RNAHP_TloopG#(0 To NumTetraLoops - 1)

'RNA enthalpy.

Private RNAHP_StackH#(0 To 3, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private RNAHP_TStackH#(0 To 3, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private RNAHP_DangleH#(0 To 1, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private RNAHP_LoopH#(0 To 2, 0 To NumZukerLoops - 1)

Private RNAHP_TloopH#(0 To NumTetraLoops - 1)

'RNA entropy.

Private RNAHP_Stacks#(0 To 3, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private RNAHP_Tstacks#(0 To 3, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private RNAHP_Dangles#(0 To 1, 0 To 3, 0 To 3, 0 To 3, 0 To 3)

Private RNAHP_Loops#(0 To 2, 0 To NumZukerLoops - 1)

Private RNAHP_Tloops#(0 To NumTetraLoops - 1)

Private Sub CalcsFromGH(Dims#, G#(), H#(), ByVal T#, S#())

'Function
'Calculate S given G, H, and T, for all elements of array S.
'Arguments
' Dims: The number of dimensions of G, H, and S.
' G: The given free energies.
' H: The given enthalpies.
' T: The temperature at which free energy was calculated.
' S: The returned entropies.

'Notes
' i. The reason for this functions existence is that Zuker provides
' G at 37, and H. S is needed to change temperatures.

Dim I#, J#, K#, L#

T = T + 273.15

Select Case Dims

Case 1

For I = LBound(G) To UBound(G)

S(I) = (G(I) - H(I)) / T

Next I

Case 2

For I = LBound(G, 1) To UBound(G, 1)

For J = LBound(G, 2) To UBound(G, 2)

S(I, J) = (G(I, J) - H(I, J)) / T

Next J

Next I

Case 3

For I = LBound(G, 1) To UBound(G, 1)

For J = LBound(G, 2) To UBound(G, 2)

For K = LBound(G, 3) To UBound(G, 3)

For L = LBound(G, 4) To UBound(G, 4)

S(I, J, K, L) = (G(I, J, K, L) - H(I, J, K, L)) / T

Next L

Next K

Next J

Next I

End Select

```

End Sub

Private Sub CalcFromHs(Dims$, H$(1), S$(1), ByVal T$, G$(1))
'Function
'Calculate G given H, S, and T, for all elements of array G.
'Arguments
' Dims: The number of dimensions of G, H, and S.
' H: The given enthalpies.
' S: The given entropies.
' T: The temperature at which free energy is required.
' S: The returned free energies.
Dim I$, J$, K$, L$
T = T + 273.15
Select Case Dims
Case 1
For I = LBound(G) To UBound(G)
G(I) = 0.1 * Cdbl(Cint(10# * (H(I) + S(I) * T)))
Next I
Case 2
For I = LBound(G, 1) To UBound(G, 1)
For J = LBound(G, 2) To UBound(G, 2)
G(I, J) = 0.1 * Cdbl(Cint(10# * (H(I, J) + S(I, J) * T)))
Next J
Next I
Case 4
For I = LBound(G, 1) To UBound(G, 1)
For J = LBound(G, 2) To UBound(G, 2)
For K = LBound(G, 3) To UBound(G, 3)
For L = LBound(G, 4) To UBound(G, 4)
G(I, J, K, L) = 0.1 * Cdbl(Cint(10# * (H(I, J, K, L) + S(I, J, K, L) * T)))
Next L
Next K
Next J
Next I
End Sub

Private Function DNA_BestHairpin$(Seq$(1))
'Function
'Compute the most stable hairpin for the given sequence.
'Arguments
' Seq: The numerical representation of the sequence.
Dim ThisHairpin$
Dim NumBases$
Dim FP$, TP$
DNA_BestHairpin = 1000#
NumBases = UBound(Seq) - LBound(Seq) + 1
For FP = 0 To NumBases - ZMinLoop - 2
For TP = FP + ZMinLoop + 1 To NumBases - 1
ThisHairpin = DNA_Hairpin(Seq, FP, TP)
If ThisHairpin < DNA_BestHairpin Then DNA_BestHairpin = ThisHairpin
Next TP
Next FP
End Function

```

```

Next TP
Next FP
End Function

Private Function RNA_BestHairpin$(Seq$(1))
'Function
'Compute the most stable hairpin for the given sequence.
'Arguments
' Seq: The numerical representation of the sequence.
Dim ThisHairpin$
Dim NumBases$
Dim FP$, TP$
RNA_BestHairpin = 1000#
NumBases = UBound(Seq) - LBound(Seq) + 1
For FP = 0 To NumBases - ZMinLoop - 2
For TP = FP + ZMinLoop + 1 To NumBases - 1
ThisHairpin = RNA_Hairpin(Seq, FP, TP)
If ThisHairpin < RNA_BestHairpin Then RNA_BestHairpin = ThisHairpin
Next TP
Next FP
End Function

Private Function RNA_Hairpin$(Seq$(1), ByVal FP$, ByVal TP$)
'Function
'Calculate the free energy of a RNA hairpin
'Arguments
' Seq: The sequence, in numeric representation.
' FP: The index of the 5'-end base that closes the loop.
' TP: The index of the 3' end base that closes the loop.
Dim Energies$(1)
Dim NumBases, StemLength, S$
'Check whether this pair can close the loop.
RNA_Hairpin = 1000#
Select Case Seq(FP)
Case 0
If Seq(TP) < 3 Then Exit Function
Case 1
If Seq(TP) < 2 Then Exit Function
Case 2
If (Seq(TP) < 1 And Seq(TP) < 3) Then Exit Function
Case 3
If (Seq(TP) < 0 And Seq(TP) < 2) Then Exit Function
End Select
'Calculate sequence size.
NumBases = UBound(Seq) - LBound(Seq) + 1
'setup to store all possible energies.

```

```

ReDim Energies(UBound(Seq) - LBound(Seq) + 1)

'Start the energy calculations with the loop.
Energies(0) = RNAHP_LoopG(2, TP - FP - 2)

'Add the stacking interaction of the first mismatch in the loop.
Energies(0) = Energies(0) + RNAHP_TstackG(Seq(FP), Seq(FP + 1), Seq(TP), Seq(TP
- 1))

'Loop over stem members
S = 0
Do While True
    FP = FP - 1
    TP = TP + 1

    'Check that there are still bases to process.
    If (FP >= 0) And (TP < NumBases) Then

        'Check that we are still in the helix.
        Select Case Seq(FP)
            Case 0
                If Seq(TP) <> 3 Then Exit Do
            Case 1
                If Seq(TP) <> 2 Then Exit Do
            Case 2
                If Seq(TP) <> 1 And Seq(TP) <> 3) Then Exit Do
            Case 3
                If Seq(TP) <> 0 And Seq(TP) <> 2) Then Exit Do
        End Select

        'Add the next stacking term
        Energies(S + 1) = Energies(S) + RNAHP_stackG(Seq(FP), Seq(FP + 1),
        Seq(TP), Seq(TP - 1))

        'Record the energy if the helix breaks here.
        Energies(S) = Energies(S) + RNAHP_DangleG(1, Seq(TP - 1), Seq(FP + 1),
        Seq(FP)) + RNAHP_DangleG(0, Seq(TP - 1), Seq(TP), Seq(FP + 1))

        S = S + 1
    Else
        Exit Do
    End If
Loop

'Add dangles, if they exist.
If (FP >= 0) Then Energies(S) = Energies(S) + RNAHP_DangleG(1, Seq(TP - 1),
Seq(FP + 1), Seq(FP))
If (TP < NumBases) Then Energies(S) = Energies(S) + RNAHP_DangleG(0, Seq(TP -
1), Seq(TP), Seq(FP + 1))

'Find the minimum energy.
RNA_Hairpin = 1000#
StemLength = S
For S = 0 To StemLength
    If Energies(S) < RNA_Hairpin Then RNA_Hairpin = Energies(S)
Next S

```

```

End Function

Private Function DNA_Hairpin$(Seq$( ), ByVal FP$, ByVal TP$)
'-----
'Function
' Calculate the free energy of a DNA hairpin
' Arguments
' Seq: The sequence, in numeric representation.
' FP: The index of the 5'-end base that closes the loop.
' TP: The index of the 3'-end base that closes the loop.
'-----
Dim Energies$( )
Dim NumBases, StemLength, S#

'Check whether this pair can close the loop.
DNA_Hairpin = 1000#
Select Case Seq(FP)
    Case 0
        If Seq(TP) <> 3 Then Exit Function
    Case 1
        If Seq(TP) <> 2 Then Exit Function
    Case 2
        If Seq(TP) <> 1 And Seq(TP) <> 3) Then Exit Function
    Case 3
        If Seq(TP) <> 0 And Seq(TP) <> 2) Then Exit Function
End Select

'Calculate sequence size.
NumBases = UBound(Seq) - LBound(Seq) + 1

'Setup to store all possible energies.
ReDim Energies(UBound(Seq) - LBound(Seq) + 1)

'Start the energy calculations with the loop.
Energies(0) = DNAHP_LoopG(2, TP - FP - 2)

'Add the stacking interaction of the first mismatch in the loop.
Energies(0) = Energies(0) + DNAHP_TstackG(Seq(FP), Seq(FP + 1), Seq(TP), Seq(TP
- 1))

'Loop over stem members
S = 0
Do While True
    FP = FP - 1
    TP = TP + 1

    'Check that there are still bases to process.
    If (FP >= 0) And (TP < NumBases) Then

        'Check that we are still in the helix.
        Select Case Seq(FP)
            Case 0
                If Seq(TP) <> 3 Then Exit Do
            Case 1
                If Seq(TP) <> 2 Then Exit Do
            Case 2
                If Seq(TP) <> 1 And Seq(TP) <> 3) Then Exit Do
        End Select

        'Add the next stacking term
        Energies(S + 1) = Energies(S) + RNAHP_stackG(Seq(FP), Seq(FP + 1),
        Seq(TP), Seq(TP - 1))

        'Record the energy if the helix breaks here.
        Energies(S) = Energies(S) + RNAHP_DangleG(1, Seq(TP - 1), Seq(FP + 1),
        Seq(FP)) + RNAHP_DangleG(0, Seq(TP - 1), Seq(TP), Seq(FP + 1))

        S = S + 1
    Else
        Exit Do
    End If
Loop

'Add dangles, if they exist.
If (FP >= 0) Then Energies(S) = Energies(S) + RNAHP_DangleG(1, Seq(TP - 1),
Seq(FP + 1), Seq(FP))
If (TP < NumBases) Then Energies(S) = Energies(S) + RNAHP_DangleG(0, Seq(TP -
1), Seq(TP), Seq(FP + 1))

'Find the minimum energy.
RNA_Hairpin = 1000#
StemLength = S
For S = 0 To StemLength
    If Energies(S) < RNA_Hairpin Then RNA_Hairpin = Energies(S)
Next S

```



```

Case 3
If (Seq(TP) <> 0 And Seq(TP) <> 2) Then Exit Do
End Select
'Add the next stacking term
Energies(s + 1) = Energies(s) + DNAHP_StackG(seq(TP), Seq(TP + 1),
Seq(TP), Seq(TP - 1))
'Record the energy if the helix breaks here.
Energies(s) = Energies(s) + DNAHP_DangleG(1, Seq(TP - 1), Seq(TP + 1),
Seq(TP)) + DNAHP_DangleG(0, Seq(TP - 1), Seq(TP), Seq(TP + 1))
S = S + 1
Else
Exit Do
End If
Loop
'Add dangles, if they exist.
If (FP >= 0) Then Energies(s) = Energies(s) + DNAHP_DangleG(1, Seq(TP - 1),
Seq(TP + 1), Seq(FP))
If (FP < Numbases) Then Energies(s) = Energies(s) + DNAHP_DangleG(0, Seq(TP -
1), Seq(TP), Seq(FP + 1))
'Find the minimum energy.
DNA_Hairpin = 1000#
StemLength = S
For S = 0 To StemLength
If Energies(s) < DNA_Hairpin Then DNA_Hairpin = Energies(s)
Next S
End Function

```

Public Sub InitThermoPars()

```

'Function:
' Initialize all the fixed parameters used for Thermo calculations.
'Notes:
' 1. These parameters are essentially constants, but are placed in
' arrays for ease-of-use. Thus this routine must be called by
' Main to initialize everything.
' 2. The nearest-neighbor parameters are accessed by indexing by
' the first then second base, with A->0, C->1, G->2, T->3
' 3. The hairpin parameters are loaded from files by this routine,
' and the entropy matrices filled in. After this operation, the
' free energy matrices may be over-written at any time.

```

```

Dim FileNames
Dim fooz()
'files opened for Zuker parameters.
'bitbucket.

```

```

'Initialize the DNA Duplex enthalpy and entropy.
'These parameters are from Santalucia et al. Biochemistry, v. 35, pp 3555.
'as found by PKW. DuplexH is in kcal/mol, DuplexS in cal/mol/deg K.

```

```

DNA_DuplexH(0, 0) = -8.4
DNA_DuplexH(0, 1) = -8.6
DNA_DuplexH(0, 2) = -6.1
DNA_DuplexH(0, 3) = -6.5
DNA_DuplexH(1, 0) = -7.4
DNA_DuplexH(1, 1) = -6.7
DNA_DuplexH(1, 2) = -10.1
DNA_DuplexH(1, 3) = DNA_DuplexH(0, 2)
DNA_DuplexH(2, 0) = -7.7
DNA_DuplexH(2, 1) = -11.1
DNA_DuplexH(2, 2) = DNA_DuplexH(1, 1)
DNA_DuplexH(2, 3) = DNA_DuplexH(0, 1)
DNA_DuplexH(3, 0) = -6.3
DNA_DuplexH(3, 1) = DNA_DuplexH(2, 0)
DNA_DuplexH(3, 2) = DNA_DuplexH(1, 0)
DNA_DuplexH(3, 3) = DNA_DuplexH(0, 0)

```

DNA_EndTAH = 0.4

```

DNA_DuplexS(0, 0) = -23.6
DNA_DuplexS(0, 1) = -23
DNA_DuplexS(0, 2) = -16.1
DNA_DuplexS(0, 3) = -18.8
DNA_DuplexS(1, 0) = -19.3
DNA_DuplexS(1, 1) = -15.6
DNA_DuplexS(1, 2) = -25.5
DNA_DuplexS(1, 3) = DNA_DuplexS(0, 2)
DNA_DuplexS(2, 0) = -20.3
DNA_DuplexS(2, 1) = -28.4
DNA_DuplexS(2, 2) = DNA_DuplexS(1, 1)
DNA_DuplexS(2, 3) = DNA_DuplexS(0, 1)
DNA_DuplexS(3, 0) = -18.5
DNA_DuplexS(3, 1) = DNA_DuplexS(2, 0)
DNA_DuplexS(3, 2) = DNA_DuplexS(1, 0)
DNA_DuplexS(3, 3) = DNA_DuplexS(0, 0)

```

```

DNA_InitGCS = -5.9
DNA_InitATS = -9
DNA_Selfs = -1.4

```

```

'Initialize the DNA/RNA duplex parameters.
'Parameters found in Sugimoto et al., Biochemistry,
'v. 34, pp. 11,211-11,216 (1995).

```

```

'Note carefully, these numbers are given for the sequence of
'the DNA strand i.e. we assume that the probe is DNA, the target
'is RNA, and the sequence of the probe is given to the routine.
'For example, the parameter in (0,0) is for dAA/rTT.

```

```

DR_DuplexH(0, 0) = -11.5
DR_DuplexH(0, 1) = -7.8
DR_DuplexH(0, 2) = -7.8
DR_DuplexH(0, 3) = -8.3
DR_DuplexH(1, 0) = -10.4
DR_DuplexH(1, 1) = -12.8
DR_DuplexH(1, 2) = -16.3
DR_DuplexH(1, 3) = -9.1
DR_DuplexH(2, 0) = -8.6

```

```

DR_DuplexH(2, 1) = -8#
DR_DuplexH(2, 2) = -9.3
DR_DuplexH(2, 3) = -5.9
DR_DuplexH(3, 0) = -7.8
DR_DuplexH(3, 1) = -5.5
DR_DuplexH(3, 2) = -9#
DR_DuplexH(3, 3) = -7.8

DR_InitH = 1.9

DR_Duplex(0, 0) = -36.4
DR_Duplex(0, 1) = -21.6
DR_Duplex(0, 2) = -19.7
DR_Duplex(0, 3) = -23.9
DR_Duplex(1, 0) = -28.4
DR_Duplex(1, 1) = -31.9
DR_Duplex(1, 2) = -47.1
DR_Duplex(1, 3) = -23.5
DR_Duplex(2, 0) = -22.9
DR_Duplex(2, 1) = -17.1
DR_Duplex(2, 2) = -23.2
DR_Duplex(2, 3) = -12.3
DR_Duplex(3, 0) = -23.2
DR_Duplex(3, 1) = -13.5
DR_Duplex(3, 2) = -26.1
DR_Duplex(3, 3) = -21.9

DR_Inits = -3.9

'Load the zucker hairpin parameters.
FileName = App.Path & "\stack.datd.pw"
ReadZucker FileName, "Stack", DNAHP_StackG, foo
FileName = App.Path & "\stack.dhd.pw"
ReadZucker FileName, "Stack", DNAHP_StackH, foo
CalcSFromGH 4, DNAHP_StackG, DNAHP_StackH, 37#, DNAHP_Stacks

FileName = App.Path & "\tstack.datd.pw"
ReadZucker FileName, "Stack", DNAHP_StackG, foo
FileName = App.Path & "\tstack.dhd.pw"
ReadZucker FileName, "Stack", DNAHP_StackH, foo
CalcSFromGH 4, DNAHP_StackG, DNAHP_StackH, 37#, DNAHP_Stacks

FileName = App.Path & "\dangle.datd.pw"
ReadZucker FileName, "Dangle", DNAHP_DangleG, foo
FileName = App.Path & "\dangle.dhd.pw"
ReadZucker FileName, "Dangle", DNAHP_DangleH, foo
CalcSFromGH 4, DNAHP_DangleG, DNAHP_DangleH, 37#, DNAHP_Dangles

FileName = App.Path & "\loop.datd.pw"
ReadZucker FileName, "Loop", DNAHP_LoopG, zLoopLengths
FileName = App.Path & "\loop.dhd.pw"
ReadZucker FileName, "Loop", DNAHP_LoopH, zLoopLengths
CalcSFromGH 2, DNAHP_LoopG, DNAHP_LoopH, 37#, DNAHP_Loops

FileName = App.Path & "\tetra.datd.pw"
ReadZucker FileName, "TetraLoop", DNAHP_TetraLoopG, zTetraLoops
FileName = App.Path & "\tetra.dhd.pw"
ReadZucker FileName, "TetraLoop", DNAHP_TetraLoopH, zTetraLoops

```

```

CalcSFromGH 1, DNAHP_TloopG, DNAHP_TloopH, 37#, DNAHP_Tloops

FileName = App.Path & "\stack.dat.pw"
ReadZucker FileName, "Stack", RNAHP_StackG, foo
FileName = App.Path & "\stack.dh.pw"
ReadZucker FileName, "Stack", RNAHP_StackH, foo
CalcSFromGH 4, RNAHP_StackG, RNAHP_StackH, 37#, RNAHP_Stacks

FileName = App.Path & "\tstack.dat.pw"
ReadZucker FileName, "Stack", RNAHP_TstackG, foo
FileName = App.Path & "\tstack.dh.pw"
ReadZucker FileName, "Stack", RNAHP_TstackH, foo
CalcSFromGH 4, RNAHP_TstackG, RNAHP_TstackH, 37#, RNAHP_Tstacks

FileName = App.Path & "\dangle.dat.pw"
ReadZucker FileName, "Dangle", RNAHP_DangleG, foo
FileName = App.Path & "\dangle.dh.pw"
ReadZucker FileName, "Dangle", RNAHP_DangleH, foo
CalcSFromGH 4, RNAHP_DangleG, RNAHP_DangleH, 37#, RNAHP_Dangles

FileName = App.Path & "\loop.dat.pw"
ReadZucker FileName, "Loop", RNAHP_LoopG, zLoopLengths
FileName = App.Path & "\loop.dh.pw"
ReadZucker FileName, "Loop", RNAHP_LoopH, zLoopLengths
CalcSFromGH 2, RNAHP_LoopG, RNAHP_LoopH, 37#, RNAHP_Loops

FileName = App.Path & "\tetra.dat.pw"
ReadZucker FileName, "TetraLoop", RNAHP_TetraLoopG, zTetraLoops
FileName = App.Path & "\tetra.dh.pw"
ReadZucker FileName, "TetraLoop", RNAHP_TetraLoopH, zTetraLoops
CalcSFromGH 1, RNAHP_TloopG, RNAHP_TloopH, 37#, RNAHP_Tloops

End Sub

Private Sub ReadZucker(FileName$, ZType$, Param$( ), FirstCol$( ))
'Function
' Read thermodynamic parameters for structure calculations.
'Arguments
' FileName: The data file.
' ZType: The type of parameter array (see below).
' Param: The parameter array to be filled in.
' FirstCol: The values from the first column.
'Notes
' 1. This routine reads data files originally designed to be read by
' MFOLD; the files must be edited to place a # sign on lines that
' do not contain data, and to replace all the " " entries by 0.
' 2. There are several types of parameter files. Type "stack"
' gives the parameters for stacking one base pair over another,
' and thus has 256 entries. Type "dangle" gives the parameters
' for dangling a base over a pair, and thus has 64 parameters.
' The other types are particular to the parameters being
' presented.
' 3. The FirstCol argument is used only by Loop and TetraLoop types.
' 4. In all these files, the order of bases is ACGU-20,1,2,3
'-----

```

```

Dim ZFile$
Dim FileLength$
Dim LineLength$
Dim WordLength$
Dim FileStr$
Dim LineStr$
Dim WordStr$
Dim A$, B$, X$, Y$, C$, T$ 'loop indices
Dim NumSeq$(0 To 3) 'numeric representation of tetraloop sequences

'Open the file, read it in.
ZFile = FreeFile
Open FileName For Binary As #ZFile
FileLength = FileLen(FileName)
FileStr = Input(FileLength, #ZFile)

'Choose how to parse the file.
Select Case ZType

```

```

Case "Stack"
'These entries represent stacking of one base pair over another:
. 5'-AX-3'.
. 3'-BY-5'.
'In each row, Y varies fast, B varies slow.
'By row, X varies fast, and A varies slowly.
'Index order used is param(A,X,B,Y).
For A = 0 To 3

```

```

For X = 0 To 3
NextLine FileStr, FileLength, LineStr, LineLength, "#"
For B = 0 To 3
For Y = 0 To 3
NextWord LineStr, LineLength, WordStr, WordLength
Param(A, X, B, Y) = Cdbl(WordStr)
Next Y
Next B
Next X
Next A

```

```

Case "Dangle"
'These entries represent danglers of 3' ends, then the 5' ends:
. 5'-AX-3'.
. 3'-B -5' is a 3' dangle
'and
. 5'-A -3'.
. 3'-BY-5' is a 5' dangle.
'In each row, X or Y varies fast, and B varies slowly.
'By row, A varies.
'Index order used is param(0,A,X,B) for 3', and param(1,A,B,Y) for the 5'
dangles.
For A = 0 To 3

```

```

NextLine FileStr, FileLength, LineStr, LineLength, "#"
For X = 0 To 3
For B = 0 To 3
NextWord LineStr, LineLength, WordStr, WordLength
Param(0, A, X, B) = Cdbl(WordStr)
Next X
Next B

```

```

Next A = 0 To 3
NextLine FileStr, FileLength, LineStr, LineLength, "#"
For B = 0 To 3
For Y = 0 To 3
NextWord LineStr, LineLength, WordStr, WordLength
Param(1, A, B, Y) = Cdbl(WordStr)
Next Y
Next B
Next A

Case "Loop"
'These entries represent internal=0, bulge=1, and hairpin=2 loops.
'By row, loop length L varies.
'Index order used is param(loop type, L)
'FirstCol is filled from the first column, representing loop lengths.
For L = 0 To NumZuckerLoops - 1
NextLine FileStr, FileLength, LineStr, LineLength, "#"
NextWord LineStr, LineLength, WordStr, WordLength
FirstCol(L) = Clng(WordStr)
For T = 0 To 2
NextWord LineStr, LineLength, WordStr, WordLength
Param(T, L) = Cdbl(WordStr)
Next T
Next L

```

```

Case "Tetraloop"
'These represent especially stable tetraloops.
'By row, sequence of the tetraloop varies.
'FirstCol is filled with a quaternary index of the sequence.
For L = 0 To NumTetraloops - 1
NextLine FileStr, FileLength, LineStr, LineLength, "#"
NextWord LineStr, LineLength, WordStr, WordLength
DNA_Str2Num WordStr, NumSeq
FirstCol(L) = 0
For C = 0 To 3
FirstCol(L) = FirstCol(L) * 4 + NumSeq(C)
Next
NextWord LineStr, LineLength, WordStr, WordLength
Param(L) = Cdbl(WordStr)
Next L

```

```

End Select
Close ZFile
End Sub

```

Function DR_CalcDeltaH\$(ByVal Seq\$)

```

'Function
. Calculate the association enthalpy for a given RNA sequence
. with its perfect W-C DNA complement.
'Arguments
. Seq: The RNA sequence.
'Returns
. Enthalpy of association, in kcal/mole.
'Notes
. 1. Standard conditions (1M NA+) is assumed.
'History

```

29-Jul-1997: From PKW's routine of the same name. PW.

Dim Length%
Dim B%
Dim NumSeq%()
'length of the sequence
'base index for traversing the sequence
'numerical representation of the sequence

'Lower case, please.
Seq = LCase(Seq)

'create the numeric representation.
Length = Len(Seq)
ReDim NumSeq(0 To Length - 1)
DNA_Str2Num Seq, NumSeq

'Initialize with initiation deltaH
DR_CalcDeltaH = DR_InitH

'Sum the nearest neighbor values along the strand.

For B = 1 To Length - 1
DR_CalcDeltaH = DR_CalcDeltaH + DR_DuplexH(NumSeq(B - 1), NumSeq(B))
Next B

'Convert to cal.
DR_CalcDeltaH = DR_CalcDeltaH * 1000#

End Function

Function DNA_CalcDeltaH(ByVal Seq%)

'Function
' Calculate the association enthalpy for a DNA sequence
' with its perfect W-C complement.

'Arguments

' Seq: The sequence.

'Returns

' Enthalpy of association, in kcal/mole.

'Notes

' 1. Standard conditions (1M NA+) is assumed.

'History

' 29-Jul-1997: From PKW's routine of the same name. PW.

Dim Length%
Dim B%
Dim NumSeq%()
'length of the sequence
'base index for traversing the sequence
'numerical representation of the sequence

'Lower case, please.
Seq = LCase(Seq)

'Create the numeric representation.

Length = Len(Seq)
ReDim NumSeq(0 To Length - 1)
DNA_Str2Num Seq, NumSeq

'Initialize with AT end correction.

If NumSeq(0) = 3 Then DNA_CalcDeltaH = DNA_CalcDeltaH + DNA_EndTAH

If NumSeq(Length - 1) = 0 Then DNA_CalcDeltaH = DNA_CalcDeltaH + DNA_EndTAH
'Sum the nearest neighbor values along the strand.

For B = 1 To Length - 1
DNA_CalcDeltaH = DNA_CalcDeltaH + DNA_DuplexH(NumSeq(B - 1), NumSeq(B))
Next B

'Convert to cal.
DNA_CalcDeltaH = DNA_CalcDeltaH * 1000#

End Function

Function DNA_CalcDeltaH(ByVal Seq%)

'Function
' Calculate the association entropy for a sequence
' with its perfect W-C complement.

'Arguments

' Seq: The sequence

'Returns

' Entropy of association, in cal/mole/deg K.

'Notes

' 1. Parameters are derived from Santalucia et al., Biochemistry,

' V. 35, pp. 3555-3562 (1996).

' 2. Standard conditions (1M NA+) is assumed.

'History

' 29-Jul-1997: From PKW's routine of the same name. PW.

Dim Length%
Dim B%
Dim NumSeq%()
'length of the sequence
'base index for traversing the sequence
'numerical representation of the sequence

'Lower case, please.
Seq = LCase(Seq)

'Create the numeric representation.

Length = Len(Seq)
ReDim NumSeq(Length - 1)
DNA_Str2Num Seq, NumSeq

'Initialize with self-symmetry correction.

If Seq = DNA_RevComp(Seq) Then DNA_CalcDeltaH = DNA_CalcDeltaH + DNA_Selfs

'add initiation term

If Instr(Seq, "c") Or Instr(Seq, "g") Then

DNA_CalcDeltaH = DNA_CalcDeltaH + DNA_InitGCS

Else

DNA_CalcDeltaH = DNA_CalcDeltaH + DNA_InitATS

End If

'Sum the nearest neighbor values along the strand

For B = 1 To Length - 1
DNA_CalcDeltaH = DNA_CalcDeltaH + DNA_DuplexS(NumSeq(B - 1), NumSeq(B))
Next B

End Function

Function DR_CalcDeltas#(ByVal Seq\$)

```

'Function
' Calculate the association entropy for a given RNA sequence
' with its perfect W-C DNA complement.
'Arguments
' Seq: The sequence
'Returns:
' Entropy of association, in cal/mole/deg K.
'Notes
' 1. Standard conditions (1M Na+) is assumed.
'History
' 29-Jul-1997: From PKW's routine of the same name. PW.

```

```

Dim Length%
Dim B%
Dim NumSeq%()
'Length of the sequence
'base index for traversing the sequence
'numerical representation of the sequence

```

```

'Lower case, please.
Seq = LCase(Seq)

```

```

'Create the numeric representation.
Length = Len(Seq)
ReDim NumSeq(0 To Length - 1)
DNA_Str2Num Seq, NumSeq

```

```

'Begin with initiation term.
DR_CalcDeltas = DR_inits

```

```

'Sum the nearest neighbor values along the strand.
For B = 1 To Length - 1
DR_CalcDeltas = DR_CalcDeltas + DR_Duplexs(NumSeq(B - 1), NumSeq(B))
Next B

```

End Function

Public Sub DNA_CalcAllTM(Seq\$(), tmp As cTMPars, TM#())

```

'Function
' Calculate the melting temperature of an array of oligos with their
' perfect W-C complements.
'Arguments
' Seq: The sequences
' TMP: An instance of the parameter class for TM calculations.
'Returns
' Melting temperature, in deg. C.
'Notes
' 1. Concentration is used as 1s, i.e. assuming the complement
' is present at much lower concentration.
'History
' 29-Jul-1997: From PKW's routine of the same name. PW.

```

On Error Goto E

```

Dim NumSeqs%
Dim S%
'number of sequences we are working with
'index

```

```

'Determine the number of probes we are calculating TM for.
NumSeqs = UBound(Seq) + 1

```

```

'Calculate melting points
For S = 0 To NumSeqs - 1
If S = Progress.StopAt = 0 Then Progress.CheckProgress S
TM(S) = DNA_CalcTM(Seq(S), tmp.Conc)

```

```

Next S
Exit Sub

```

```

E: Debug.Print "Error in DNA_CalcAllTM"
Err.Raise Err.Number, , Err.Description
End Sub

```

Public Sub DNA_CalcClamp(Seq\$(), TClamp As cClampPars, Clamp#())

```

'Function
' Calculate the melting temperature of the clamp of a probe
' to its perfect W-C complement.
'Arguments
' Seq: The sequences
' TClamp: An instance of the parameter class for Clamp calculations.
'Returns
' Melting temperature of tightest clamp, in deg. C.

```

```

'Notes
' 1. Concentration is used as 1s, i.e. assuming the complement
' is present at much lower concentration.

```

On Error Goto E

```

Dim NumSeqs%
Dim S%, SS%
Dim SubSeq$
Dim BestTM#, ThisTM#
'number of sequences we are working with
'indices
'subsequence
'current most stable clamp

```

```

'Determine the number of probes we are calculating clamp for.
NumSeqs = UBound(Seq) + 1

```

```

'Calculate melting points.
For S = 0 To NumSeqs - 1
If S = Progress.StopAt = 0 Then Progress.CheckProgress S
BestTM = 0

```

```

For SS = TClamp.Fivep + 1 To Len(Seq(S)) - TClamp.Threep - TClamp.Length + 1
SubSeq = Mid(Seq(S), SS, TClamp.Length)
ThisTM = DNA_CalcTM(SubSeq, TClamp.Conc) + 273.15
If ThisTM > BestTM Then BestTM = ThisTM

```

```

Next SS
Clamp(S) = BestTM

```

```

Next S
Exit Sub
E: Debug.Print "Error in DNA_CalcClamp"
Err.Raise Err.Number, , Err.Description
End Sub

```

```

Public Sub DR_CalcClamp(Seqs(), TClamp As cClampParams, Clamp#())
'Function
' Calculate the melting temperature of the Clamp of a probe
' to its perfect RNA W-C complement.
'Arguments
' Seq: The sequences
' TClamp: An instance of the parameter class for Clamp calculations.
'Returns
' Melting temperature of tightest clamp, in deg. C.
'Notes
' 1. Concentration is used as 1s, i.e. assuming the complement
' is present at much lower concentration.
On Error GoTo E
Dim NumSeqs#
Dim S#, SS#
Dim SubSeq#
Dim BestTM#, ThisTM#
'number of sequences we are working with
'indices
'subsequence
'current most stable clamp
'Determine the number of probes we are calculating Clamp for.
NumSeqs = UBound(Seq) + 1
'Calculate melting points.
For S = 0 To NumSeqs - 1
    BestTM = 0
    If S - Progress.StopAt = 0 Then Progress.CheckProgress S
    For SS = TClamp.Fivep + 1 To Len(Seq(S)) - TClamp.Threep - TClamp.Length + 1
        SubSeq = Mid(Seq, SS, TClamp.Length)
        ThisTM = DR_CalcTM(subseq, TClamp.Conc) + 273.15
        If ThisTM > BestTM Then BestTM = ThisTM
    Next SS
    Clamp(S) = BestTM
Next S
Exit Sub
E: Debug.Print "Error in DR_CalcClamp"
Err.Raise Err.Number, , Err.Description
End Sub

Public Sub DNA_CalcCdGH(Seqs(), dGHP As cDGHParams, dGH#())
'Function
' Calculate the hairpin dGs of an array of oligos.
'Arguments
' Seq: The sequences.
' dGHP: An instance of the parameter class for dGH calculations.
' dGH: The hairpin dGs.
On Error GoTo E
Dim NumSeq#()
Dim S#
'numeric representation of the sequence
'index
'Recalculate all G parameter matrices at current temperature.

```

```

'Calculate dGs.
For S = 0 To UBound(seq)
    If S - Progress.StopAt = 0 Then Progress.CheckProgress S
    ReDim NumSeq(0 To Len(seq(S)) - 1)
    DNA.Str2Num Seq(S), NumSeq
    dGH(S) = DNA_BestHairpin(NumSeq)
Next S
Exit Sub
E: Debug.Print "Error in DNA_CalcGH"
Err.Raise Err.Number, , Err.Description
End Sub

Public Sub DNA_CalcdGD(seq#(), dGDP As cGDParams, dGD#())
'Function
' Calculate the duplex dGs of an array of oligos.
'Arguments
' Seq: The sequences.
' dGDP: An instance of the parameter class for dGD calculations.
' dGD: The hairpin dGs.
On Error GoTo E
Dim NumSeq#()
Dim S#
'numeric representation of the sequence
'index
'Calculate dGs.
For S = 0 To UBound(Seq)
    If S - Progress.StopAt = 0 Then Progress.CheckProgress S
    dGD(S) = (DNA_CalcDeltaH(seq(S)) - (dGDP.T + 273.15) *
    DNA_CalcDeltaS(seq(S))) / 1000#
Next S
Exit Sub
E: Debug.Print "Error in DNA_CalcGD"
Err.Raise Err.Number, , Err.Description
End Sub

Public Sub DR_CalcdGD(Seqs(), dGDP As cGDParams, dGD#())
'Function
' Calculate the duplex dGs of an array of oligos.
'Arguments
' Seq: The sequences.
' dGDP: An instance of the parameter class for dGD calculations.
' dGD: The hairpin dGs.
On Error GoTo E
Dim NumSeq#()
Dim S#
'numeric representation of the sequence
'index
'Calculate dGs.

```

```

For S = 0 To UBound(Seq)
    If S - Progress.StopAt = 0 Then Progress.CheckProgress S
    dGH(S) = (DR_CalcDeltaH(Seq(S)) - (dGDP.T + 273.15) * DR_CalcDeltaS(Seq(S)))
    / 1000#
Next S
Exit Sub
E: Debug.Print "Error in DNA.CalcGD"
Err.Raise Err.Number, , Err.Description
End Sub

Public Sub DNA_CalcDGM(Seq$(), dGMP As cdGMPars, dGM#())
    'Function
    'Calculate the MFold dGs of an array of oligos.
    'Arguments
    ' Seq: The sequences.
    ' dGMP: An instance of the parameter class for dGM calculations.
    ' dGM: The hairpin dGs.
    '-----
    On Error Goto E

Dim S# 'Index

'Calculate dGs.
For S = 0 To UBound(Seq)
    If S - Progress.StopAt = 0 Then Progress.CheckProgress S
    dGM(S) = frmMain.Mfoldx.mfold(Seq(S), Val(dGMP.T), True)
Next S
Exit Sub
E: Debug.Print "Error in DNA.CalcDGM"
Err.Raise Err.Number, , Err.Description
End Sub

Public Sub RNA_CalcDGH(Seq$(), dGHP As cdGHPars, dGH#())
    'Function
    'Calculate the hairpin dGs of an array of oligos.
    'Arguments
    ' Seq: The sequences.
    ' dGHP: An instance of the parameter class for dGH calculations.
    ' dGH: The hairpin dGs.
    '-----
    On Error Goto E

Dim NumSeq#() 'numeric representation of the sequence
Dim S# 'Index

'Recalculate all G parameter matrices at current temperature.
CalcGFromHS 4, RNAHP_StackH, RNAHP_Stacks, dGHP.T, RNAHP_StackG
CalcGFromHS 4, RNAHP_TStackH, RNAHP_TStacks, dGHP.T, RNAHP_TStackG
CalcGFromHS 4, RNAHP_DangleH, RNAHP_Dangles, dGHP.T, RNAHP_DangleG
CalcGFromHS 2, RNAHP_LoopH, RNAHP_Loops, dGHP.T, RNAHP_LoopG
CalcGFromHS 1, RNAHP_TloopH, RNAHP_Tloops, dGHP.T, RNAHP_TloopG

'Calculate dGs.

```

```

For S = 0 To UBound(Seq)
    If S - Progress.StopAt = 0 Then Progress.CheckProgress S
    ReDim NumSeq(0 To Len(Seq(S)) - 1)
    DNA_Str2Num Seq(S), NumSeq
    dGH(S) = RNA_BestHairpin(NumSeq)
Next S
Exit Sub
E: Debug.Print "Error in RNA.CalcDGH"
Err.Raise Err.Number, , Err.Description
End Sub

Public Sub RNA_CalcDGM(Seq$(), dGMP As cdGMPars, dGM#())
    'Function
    'Calculate the MFold dGs of an array of oligos.
    'Arguments
    ' Seq: The sequences.
    ' dGMP: An instance of the parameter class for dGM calculations.
    ' dGM: The MFold dGs.
    '-----
    On Error Goto E

Dim S# 'Index

'Calculate dGs.
For S = 0 To UBound(Seq)
    If S - Progress.StopAt = 0 Then Progress.CheckProgress S
    dGM(S) = frmMain.Mfoldx.mfold(Seq(S), Val(dGMP.T), False)
Next S
Exit Sub
E: Debug.Print "Error in RNA.CalcGM"
Err.Raise Err.Number, , Err.Description
End Sub

Public Sub DR_CalcAllTM(Seq$(), tmp As cTMPars, TM#())
    'Function
    'Calculate the melting temperature of an array of DNA probes with their
    ' perfect W-C RNA target complements.
    'Arguments
    ' Seq: The sequences.
    ' TMP: An instance of the parameter class for TM calculations.
    'Returns
    ' Melting temperature, in deg. C.
    'Notes
    ' 1. Concentration is used as is, i.e. assuming the complement
    ' is present at much lower concentration.
    'History
    ' 29-Jul-1997: From PKW's routine of the same name. PW.
    '-----
    On Error Goto E

Dim NumSeq# 'number of sequences we are working with
Dim S# 'index

'Determine the number of probes we are calculating TM for.
NumSeqs = UBound(Seq) - LBound(Seq) + 1

```

```
'Calculate melting points
For S = 0 To NumSeqs - 1
    If S = Progress.StopAt = 0 Then Progress.CheckProgress S
    TM(S) = DR_CalcTM(Seq(S), tmp.Conc)
Next S
Exit Sub
E: Debug.Print "Error in DR_CalcAllTM"
Err.Raise Err.Number, , Err.Description
End Sub

Public Function DR_CalcTM#(Seq$, Conc#)
    'Function
    'Calculate the TM of one DNA/RNA duplex. Sequence given is DNA.
    DR_CalcTM = DR_CalcDeltaH(Seq) / (DR_CalcDeltaS(Seq) + RGas * Log(Conc)) -
    273.15
End Function
```

```
Public Function DNA_CalcTM#(Seq$, Conc#)
    'Function
    'Calculate the TM of one DNA/DNA duplex.
    DNA_CalcTM = DNA_CalcDeltaH(Seq) / (DNA_CalcDeltaS(Seq) + RGas * Log(Conc)) -
    273.15
End Function
```

```
Attribute VB_Name = "Utilities"
Option Explicit
```

```
'A quicksort routine for doubles, based on NR code.
Private Declare Function vbSort2 Lib "vb5nrdll.dll" _
    (ByVal N$, ByVal Vector$, ByVal Index$) As Long
```

```
Public Sub CalcRun(Pos$, RP As cRunParms, Run$())
```

```
    'Function
    'Calculate the runs in a set of positions.
```

```
    'Arguments
    ' Pos: The positions.
    ' RP: The run parameters.
    ' Run: The returned run.
```

```
    'Method
    'Begin by looping over all positions, extending down over all consecutive
    'entries, then up over all consecutive entries. Then mark in the active
    'position the length of run found. Next, eliminate all runs that are
    'too short. Finally, for each run, mark the requested number of members
    'of the run, at the requested spacing, as being members. If the run is
    'shorter than can accommodate the requested number of members, fill in as
    'many as possible, while still preserving the requested spacing.
```

```
    'Notes
    ' 1. Run must be allocated and sized correctly by the caller.
```

```
    Dim P$
    Dim EndUp$
    Dim EndDown$
    'index
    'end of the run, counting up
    'end of the run, counting down
```

```
Dim StepDir$
Dim M$
'current step direction
'index

'Loop over positions.
For P = 0 To UBound(Pos)
    EndDown = P
    Do While ((P - EndDown) = (Pos(P) - Pos(EndDown)))
        EndDown = EndDown - 1
        If EndDown = -1 Then Exit Do
    Loop
    EndDown = EndDown + 1
    EndUp = P
    Do While ((EndUp - P) = (Pos(EndUp) - Pos(P)))
        EndUp = EndUp + 1
        If EndUp = UBound(Pos) + 1 Then Exit Do
    Loop
    EndUp = EndUp - 1
    Run(P) = EndUp - EndDown + 1
Next P

'Remove runs that are too short.
For P = 0 To UBound(Run)
    If Run(P) < RP.Min Then Run(P) = 0
Next P

'Pick out requested elements from each run.
P = 0
Do While P < UBound(Run)
    'Find next run.
    Do While Run(P) = 0 And P < UBound(Run)
        P = P + 1
    Loop
    'Record the ends, step to the middle (to the right of middle for even
    'lengths).
    EndDown = P
    EndUp = P + Run(P) - 1
    P = Fix((EndDown + EndUp + 1) / 2)
    'Mark elements, stepping down from the middle first. This ensures that all
    'elements get marked for even lengths, spacing=1.
    StepDir = -1
    For M = 1 To RP.Num
        Run(P) = -Run(P)
        P = P + StepDir * M * RP.Spacing
        If (P < EndDown Or P > EndUp) Then Exit For
        StepDir = StepDir * -1
    Next M
    'Move over this run.
    P = EndUp + 1
Loop

'Convert marked elements back to run length.
For P = 0 To UBound(Run)
    If Run(P) > 0 Then Run(P) = 0
    If Run(P) < 0 Then Run(P) = -Run(P)
```


Next P

End Sub

Public Function IsGoodRS(RS As Recordset) As Boolean

```

Function
    Check whether it is possible to access the records of
    recordset connected to the sequence, proberset, or probes table.
    No current record need exist for this function to return True.

```

```

On Error GoTo Err
IsGoodRS = False
If IsNull(RS) Then Exit Function
If (RS.EOF = True And RS.BOF = True) Then Exit Function
IsGoodRS = True
Err:
End Function

```

Public Function NumRecords(RS As Recordset)

```

Function
    Count the number of records in the recordset.
Notes
    A side effect of this routine is to position the recordset at the
    first record.

```

```

RS.MoveLast
RS.MoveFirst
NumRecords = RS.RecordCount
End Function

```

Public Sub UnpackSequence(ByRef SeqStr As String, ByRef NumStr As String)

```

SeqStr is assumed to hold a packed sequence.
Dim NewSeqStr$
Dim Bases$, I$
NumStr = ""

```

```

Bases = 1
Do While Len(SeqStr) - Bases > 60
    NumStr = NumStr & Format(Format(Bases, "#####"), "#####") & vbCrLf
    For I = Bases To Bases + 59
        NewSeqStr = NewSeqStr & Mid(SeqStr, I, 1)
        If (I Mod 10 = 0) Then NewSeqStr = NewSeqStr & " "
    Next I
    Bases = Bases + 60
    NewSeqStr = NewSeqStr + vbCrLf

```

```

Loop
If Bases <> Len(SeqStr) Then
    NumStr = NumStr & Format(Format(Bases, "#####"), "#####")
    For I = Bases To Len(SeqStr)
        NewSeqStr = NewSeqStr & Mid(SeqStr, I, 1)
        If (I Mod 10 = 0) Then NewSeqStr = NewSeqStr & " "
    Next I

```

```

End If
SeqStr = NewSeqStr
End Sub

```

```

Public Function StrField(str$, lineLen$, ByVal FieldLen$, Optional Extra$)
'Returns a leading field of a string, and shortens the string
'by the fieldlength + extra.
'If the line is too short, as much of the requested field as possible
'will be returned.

```

```

If FieldLen > lineLen Then FieldLen = lineLen
StrField = TrimLeft(str, FieldLen)
lineLen = lineLen - FieldLen
If Not IsMissing(Extra) Then
    lineLen = lineLen - Extra
    If lineLen < 0 Then lineLen = 0

```

```

End If
str = Right(str, lineLen)
End Function

```

Public Sub NextWord(str\$, StrLen\$, Word\$, WordLen\$)

```

'Strips leading spaces, copies leading word, reduces line length
Do While (Mid(str, 1, 1) = " ") And (StrLen > 0)
    StrLen = StrLen - 1
    str = Right(str, StrLen)

```

```

Loop
If Instr(str, " ") <> 0 Then
    WordLen = Instr(str, " ") - 1
    Else
        WordLen = StrLen

```

```

End If
If WordLen <> 0 Then
    Word = Mid(str, 1, WordLen)
    StrLen = StrLen - WordLen
    str = Right(str, StrLen)

```

```

Else
    Word = ""
    str = ""
    StrLen = 0
    WordLen = 0

```

```

End If
End Sub

```

Public Sub NextLine(FileStr\$, FileLen\$, LineStr\$, lineLen\$, CommentChar\$)

```

'Return the next line that doesn't begin with CommentChar
lineLen = 0
Do While (lineLen = 0 And FileLen > 0)
    lineLen = Instr(FileStr, vbCrLf) - 1
    LineStr = Left(FileStr, lineLen)
    FileLen = FileLen - lineLen - 2

```

```

FileStr = Right(FileStr, FileLen)
If (Mid$(LineStr, 1, 1) = CommentChar) Then lineLen = 0
Loop

```

Public Sub QuickSort(X\$(), Index\$(), LB\$, UB\$)

```

'QuickSorts the array X into ascending order,
'simultaneously sorting I to provide a sort index.
Dim N$, foot
N = UB - LB + 1
foo = vb5sort2(N, X(LB), Index(LB))

```

```

End Sub

Public Function BinarySearch$(Vector$( ), Value#, U$, L$)
'-----
'Function
' Find the index of a value in a sorted array.
'-----
Do While U <> L
    BinarySearch = (U + L) / 2
    If Vector(BinarySearch) = Value Then Exit Function
    If Vector(BinarySearch) > Value Then
        U = BinarySearch
    Else
        L = BinarySearch
    End If
Loop
End Function

Public Function PackSequences$(SeqStr$)
Dim NewStr$, Ch$
Dim I%
For I = 1 To Len(SeqStr)
    Ch = Mid$(SeqStr, I, 1)
    If Instr("ACGTUacgtu", Ch) Then PackSequence = PackSequence & Ch
Next I
End Function

```